

### Seat to Floor Height Ranges

		Real Wheels	Casters				
			3"x1"	4" x 1"	5.5" x 1.5"	6" x 2"	8" x 2"
Rear Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"	
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"	
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	
Reverse Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"	
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"	
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	

**We are a custom company  
If you don't see it  
Please ask**