

Account # _____
 Phone _____
 Bill to _____

 Purchase Order # _____
 Account Contact _____
 Email Address _____

 Authorized By _____
 Fax # _____
 Ship To _____

 Date _____
 Client Name _____

STANDARD FEATURES

- All Walkers
 - Heavy duty frame and parts
 - Non-Tipping Base made of durable steel
 - Height Adjustable
 - Ring Style Padding Allows for Comfortable Seating
 - Lifetime warranty on base
 - 5 year warranty on all Gunnell designed products
- Infant Walker
 1. Patient Height, Approx: 27" – 36"
 2. Balance Ring Height:
 - A. Minimum: 18"
 - B. Maximum: 24"
 3. Ring Size Opening: 8 ½" x 10"
 4. Base Width: 19"
 5. Base Length: 20"
 6. Weight (lbs): 15
- Youth Walker
 1. Patient Height, Approx: 34" – 54"
 2. Balance Ring Height:
 - A. Minimum: 22"
 - B. Maximum: 35"
 3. Ring Size Opening: 12" x 10 ½"
 4. Base Width: 31"
 5. Base Length: 34"
 6. Weight (lbs): 20
- Adolescent Walker
 1. Patient Height, Approx: 36" – 56"
 2. Balance Ring Height:
 - A. Minimum: 27"
 - B. Maximum: 42"
 3. Ring Diameter: 15"
 4. Base Width: 33"
 5. Base Length: 40"
 6. Weight (lbs): 30
- Adult Walker
 1. Patient Height, Approx: 48" – 68"
 2. Balance Ring Height:
 - A. Minimum: 43"
 - B. Maximum: 58"
 3. Ring Diameter: 16"
 4. Base Width: 34"
 5. Base Length: 47"
 6. Weight (lbs): 35

✓ SUPPORT WALKERS

- | | | |
|-------------------------------|------------------------|------------|
| <input type="checkbox"/> IW-1 | Infant Walker..... | \$559.00 |
| <input type="checkbox"/> YW-1 | Youth Walker..... | \$749.00 |
| <input type="checkbox"/> AD-1 | Adolescent Walker..... | \$982.00 |
| <input type="checkbox"/> AW-1 | Adult Walker..... | \$1,199.00 |
| <input type="checkbox"/> CW-1 | Custom Walker..... | Quote |

✓ COLOR N/C

- | | | |
|---|---|---|
| <input type="checkbox"/> Pink | <input type="checkbox"/> Dark Maroon | <input type="checkbox"/> Purple Metallic |
| <input type="checkbox"/> Flat Black | <input type="checkbox"/> Sparkle Orange | <input type="checkbox"/> Green Apple Metallic |
| <input type="checkbox"/> Glitter Silver | <input type="checkbox"/> Glitter Blue | <input type="checkbox"/> Candy Blue |
| <input type="checkbox"/> Candy Purple | <input type="checkbox"/> Candy Red | <input type="checkbox"/> Sparkle Green |
| <input type="checkbox"/> Bright Yellow | <input type="checkbox"/> Blood Red | <input type="checkbox"/> Medium Blue |
| <input type="checkbox"/> ¹ Silver Vein | <input type="checkbox"/> Gloss Black | <input type="checkbox"/> Black Stardust |


✓ UPHOLSTERY COLORS

- | | | |
|--------------------------------------|---|---------|
| <input type="checkbox"/> Black | <input type="checkbox"/> ¹ Blue..... | N/C |
| <input type="checkbox"/> Custom..... | | \$49.00 |

INFANT WALKER

E0140
YOUTH WALKER

E0140
ADOLESCENT WALKER

E0140
ADULT WALKER

E0140
REPLACEMENT PARTS

- | | | |
|----------------------------------|--|----------|
| <input type="checkbox"/> YW-SB | Scissor Board for Youth Walker..... | \$99.00 |
| <input type="checkbox"/> AD-SB | Scissor Board for Adolescent Walker..... | \$159.00 |
| <input type="checkbox"/> AW-SB | Scissor Board for Adult Walker..... | \$199.00 |
| <input type="checkbox"/> S-040 | Infant Ring Pads..... | \$99.00 |
| <input type="checkbox"/> S-039 | Youth Ring Pads..... | \$109.00 |
| <input type="checkbox"/> S-087 | Adolescent Ring Pads..... | \$119.00 |
| <input type="checkbox"/> S-038-F | Adult Ring Pad Rear..... | \$59.00 |
| <input type="checkbox"/> S-038-B | Adult Ring Pad Front..... | \$59.00 |
| <input type="checkbox"/> S-042 | Infant Sling..... | \$79.00 |
| <input type="checkbox"/> S-051 | Youth Sling..... | \$89.00 |
| <input type="checkbox"/> S-088 | Adolescent Sling..... | \$99.00 |
| <input type="checkbox"/> S-041 | Adult Sling..... | \$109.00 |

RUBATEX, FLEECE OR DARLEX COVERS

We can cover any pad with Rubatex, Fleece or Darlex
Please call for quotes

NOTE: HCPCS Codes provided were determined by SADMERC and/or interpretations of the codes. Proper use of codes does not ensure coverage or payment. Please contact a Medicare professional for further assistance.

FOOTNOTES

1: Gunnell will select factory settings if no setting is made