

# Rehab MAC Order Form



Account # \_\_\_\_\_  
 Client Name \_\_\_\_\_  
 Phone \_\_\_\_\_  
 Bill to \_\_\_\_\_  
 Purchase Order # \_\_\_\_\_

Fax # \_\_\_\_\_  
 Ship To \_\_\_\_\_  
 \_\_\_\_\_  
 Date \_\_\_\_\_

## STANDARD FEATURES

- 50° Quick Release Tilt, 3° Anterior and 47° Posterior with adjustable pivot point
- 75° Seat to Back Angle Adjustable
- 3 Seat to Floor Heights to Choose from
- Wheel Locks with adjustable Tension
- Tip Leverage Bar
- Lifetime warranty on base and 1 year no question asked warranty on all Gunnell designed products.

**BASE PRICE** ....Adjustable 15 3/4" – 19 3/4" .....\$2,945.00

## FACTORY SETTINGS

√ **STANDARD COLOR** ..... N/C

Gloss Black  Green Apple Metallic  Red  Medium Blue

√ **CUSTOM COLOR** ..... N/C

Pink  Dark Maroon  Black Stardust  Sparkle Orange  
 Glitter Silver  Glitter Blue  Candy Blue  Candy Purple  
 Candy Red  Sparkle Green  Bright Yellow  Silver Vein

√ **FRAME WIDTH** ..... N/C

14 3/4"  15 3/4"  16 3/4"  17 3/4"  18 3/4"  19 3/4"

√ **SEAT DEPTH** ..... N/C

15"  16"  17"  18"  19"

√ **SEAT TO FLOOR HEIGHT** ..... N/C

<sup>1</sup> 19" Regular  17" Hemi <sup>2</sup> 15" Super-Hemi

√ **BACK CANE HEIGHT** ..... N/C

(Cane Height is measured from the Seat Rail)

<sup>1</sup> 21"  22"  23"  24"  25"  26"  27"

√ **ARMREST HEIGHT OF SEAT RAIL** ..... N/C

<sup>1</sup> 8"-12"  8"-16"

√ **FOOTREST LEGDROP** ..... N/C

(Keep Approx. Thickness of Seat Cushion in Mind When Ordering)

<sup>1</sup> 11"-17"  16"-22"

### ADJUSTABLE SEAT TO BACK ANGLE

ST Standard Tilt in Space 47° posterior to 3° anterior.... N/C

<sup>1</sup> 505 Standard 75° Adjustable seat to Back Angle  
 (85° to 160°)..... N/C

504 Extended Seat to Back Angle  
 (Opens angle up to 170°).....\$50.00

503 Quick Release Recline  
 (0° to 135° & 55° to 170°).....\$599.00

### PUSH HANDLE UPRIGHTS

<sup>1</sup> 570 Regular Height (38 1/2")..... N/C

575 Push Handle, Angle & Height Adjustable.....\$245.00

579 Height Adjustable (7").....\$139.00

## ARMREST OPTIONS

<sup>1</sup> 720S Heavy Duty, quick Release, Height & Lateral  
 Adjustable Seat Rails with Pads.....\$225.00

727X Desk Arm Length Pads.....\$225.00

<sup>6</sup> 735X Swing-Up & Height Adjustable Back  
 Canes with Pads.....\$225.00

745X Angle Adjustable Brackets(735X only).....\$129.00

728X Durable Plastic Pads.....\$39.00

441 Gel Pads.....\$120.00

## FOOTREST OPTIONS

460T 60° Swing Away Lift Off Leg rests.....\$284.75/Pair

<sup>1</sup>470T 70° Swing Away Lift Off Leg rests.....\$284.75/Pair

490T 90° Swing Away Lift Off Leg rests.....\$284.75/Pair

475T Low-shear elevating, Swing Away Lift Off  
 Leg rests With Calf Pads, pair.....\$295.00

496T Ankle adjustable footplates.....\$100.00/each

463S One Non-Elev. 90° Flip-Up Footrest  
 with Calf Pad.....\$200.00

439 One Foot Plate Pad.....\$89.00

460X Knee Angle Adjustable Hardware,  
 manual lock (for 463S only).....\$131.00

468X Ankle Angle Adjustable Hardware  
 (for 463S Only).....\$99.00

ABS-FTR Foot box (No pad).....\$195.00

432 Padded foot box.....\$295.00

428 Reinforced Footrest Bracket.....\$66.00

446 One Footrest Bumper Guard.....\$55.00

## CASTER OPTIONS

6X2SPN 6" X 2" Semi-pneumatic Casters.....\$75.00

<sup>1</sup> 556S 8" X 2" Semi-pneumatic Casters.....\$75.00

555X 5-1/2" Semi-pneumatic Casters..... N/C

## WHEEL OPTIONS

553UM 24" x 1 3/8" Urethane Tire/Composite Wheels  
 with Hand rims.....\$352.90/Pair

553M 24" x 1 3/8" Pneumatic Tire/Composite Wheels with  
 Hand rims.....\$352.90/Pair

553FM Flat Free Inserts for 553M.....\$125.00

551UM 20" x 1 3/8" Urethane Tire/Composite Wheels  
 with Hand rims..... N/C

551M 20" x 1 3/8" Pneumatic Tire/Composite Wheels with  
 Hand rims.....\$352.90/Pai

551FM Flat Free Inserts for 551M.....\$125.00

## FRAME OPTIONS

RR-WL Wheel Locks..... N/C

574 Wheel Lock Extension Handles.....\$42.00/each

569 Anti-Tips (Height & Depth Adjustable).....\$195.00

551A Quick Release Axles.....\$75.00

590 One Arm Drive  Left  Right .....\$2,500.00

TP Transportation Package.....\$250.00

## POWER OPTIONS

600 Power Actuation.....QUOTE

Tilt in Space  Lateral Tilt  Recline

591 Battery Box or Vent Tray.....\$299.00

(Specify Dimensions) Length \_\_\_\_\_" Width \_\_\_\_\_" Height \_\_\_\_\_"

\*See Footnotes at bottom of page 2

50 Enterprise Drive, Vassar, MI 48768 P: (800) 551-0055 F: (800) 794-5483 www.Gunnell-Inc.com


F010MAC-1304

# Rehab MAC Order Form


## SEAT BASES (SPECIFY DIMENSIONS)

All seats include standard attaching hardware & 1" Sunmate foam with 1/2" poly foam.


- 020S Solid Seat.....\$400.00




- 030S Anti-Thrust Seat.....\$450.00




- 040S Wedge Seat.....\$450.00



- 050S I-Modification Seat.....\$500.00



- 060S Contoured Seat.....\$525.00



- AQRSH Quick Release Hardware upcharge.....add \$25.00

- 328 2" Sunmate Foam Only.....add \$149.00

## ATTACHMENT HARDWARE


- SAH Standard Attaching Hardware.....\$140.00
- QRSB Quick Release Seat Hardware.....\$165.00

## SEAT ACCESSORIES

- 652 2" Hip Belt.....\$85.00
- 650 1.5" Hip Belt.....\$60.00
- 337 Hip Belt Pad.....\$29.00
- 634 Reinforced Seat Bracket.....\$190.00

## HIP PADS


- 656X Hip Pads (each).....\$112.00




## BACK BASES (SPECIFY DIMENSIONS)

All backs include standard attaching hardware & 1" Sunmate foam with 1/2" poly foam.


- 020B Solid Back.....\$400.00




- 030B Bi-Angular Back.....\$450.00



- 040B Curved Back.....\$450.00



- 050B I-Modification Back.....\$500.00



- AQRBH Quick Release Hardware Upcharge.....add \$25.00
- 326X 2" Sunmate Foam Only.....add \$149.00

## ATTACHMENT HARDWARE

- BAH Standard Back Hardware.....\$140.00
- QRBH Quick Release Back Hardware.....\$165.00

## SEAT ACCESSORIES

- 334 Chest Strap Ring Style.....\$97.00
- 335 Butterfly Harness.....\$137.00

## TORSO PADS

- 306 Torso Pads (pair).....\$200.00 (SPECIFY INSET & SELECT SIZE)
  - Small 2 1/2" x 4 1/4"
  - Medium 2 3/4" x 5 1/4"
  - Large 3 3/4" x 6 1/2"

### Inset

- 0.00"
- 1.00"
- 1.50"
- 2.00"

- 307 Swing Away Torso Pads (pair).....\$265.00 (SPECIFY INSET & SELECT SIZE)
  - Small 2 1/2" x 4 1/4"
  - Medium 2 3/4" x 5 1/4"
  - Large 3 3/4" x 6 1/2"


### Inset

- 0.00"
- 1.00"
- 1.50"
- 2.00"

## OCCIPITAL HEAD CRADLES (DIMENSIONS SHOWN)

All Occipital Head Cradles include attaching hardware.

- 230C Contoured Head & Neck Support.....\$295.00



- 230S Small Occipital Cradle.....\$375.00
- 231S Small Occipital Cradle w/Mesh Insert.....\$425.00



- 230M Medium Occipital Cradle.....\$395.00
- 231M Medium Occipital Cradle w/Mesh.....\$445.00



- 230L Large Occipital Cradle.....\$450.00
- 231L Large Occipital Cradle w/Mesh Insert.....\$500.00



## ABDUCTORS

- 402 Regular Abduction System.....\$225.00
- 403 Flip Down Abduction System.....\$405.00

## ADDUCTORS

- 405R Right Side Extended Hip Side.....\$65.00

Length\_\_\_\_\_

- 405L Left Side Extended Hip Side.....\$65.00

Length\_\_\_\_\_

## UPPER EXTREMITY TRAYS

- 531T Clear Tray w/Trough Style Hardware.....\$318.00
- 533T ABS Tray w/Trough Style Hardware.....\$260.00
- 534T Padded ABS Tray w/Trough Style Hardware.....\$295.00

**WE ARE A CUSTOM COMPANY.  
IF YOU DON'T SEE IT,  
PLEASE ASK!**

## FOOTNOTES

- 1: Accepts seals up to 2 1/2" wider (O.D.) than Seat rail width.
- 2: Does NOT accept seats wider than seat rail width
- 3: Comes with Anti-Tips, Handrims and Rear Wheel Axles mounted 3" forward
- 4: For 030S Seat or 668 option only.
- 5: Only available with the 030S Starboard Seat
- 6: Only available with the 120S Flat Back NO wider than the seat rail width
- 7: NOT available with the 463S footrest.
- 8: Only available with the foot box style footrest.
- 9: Measurements given are the inside dimension; add 1" for outside dimension.