



Rehab Bariatric TNT

Adjustable Bariatric Tilt-in-Space Wheelchair

Date _____ Account # _____
 Company _____
 Account Contact _____
 Bill to Address _____
 City _____ ST _____ Zip _____
 Tel: _____ Fax _____
 Email: _____

PO # _____
 Ship to _____
 Ship to Attention _____
 Ship to Address _____
 City _____ ST _____ Zip _____
 Tel: _____ Fax _____
 Comments _____

STANDARD FEATURES

- Heavy duty frame and parts (750 lbs Weight Capacity)
- Infinite adjustment 55° tilt range (default setting -10° to 45°)
- 90° - 130° of adjustable static recline in 10° increments
- 15"-19" seat depth adjustment in 1" increments
- 16"-19" frame width built in adjustment in 1/2" increments
- 15"-19" seat to floor height adjustment in 1" increments
- 17"-27" push handle height adjustment (maintain floor to handle height)
- 8" of forward to wheel base adjustment
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding

E1161 **BASE PRICE**\$2,945.00
 K0108 **BT-HDP EXTRA HEAVY DUTY PACKAGE**\$800.00

FRAME OPTIONS

CLASSIC COLOR N/C

- | | | | |
|--|--|--|--|
| <input type="checkbox"/> Gloss Black | <input type="checkbox"/> Graphite Black | <input type="checkbox"/> Black Wrinkle | <input type="checkbox"/> Lime Green |
| <input type="checkbox"/> Forest Green | <input type="checkbox"/> Cornflower Blue | <input type="checkbox"/> Baby Blue | <input type="checkbox"/> Midnight Blue |
| <input type="checkbox"/> Tangerine | <input type="checkbox"/> Ruby Red | <input type="checkbox"/> Traffic Red | <input type="checkbox"/> Maroon |
| <input type="checkbox"/> Sunny Yellow | <input type="checkbox"/> Bubblegum Pink | <input type="checkbox"/> Raspberry | <input type="checkbox"/> Blue Violet |
| <input type="checkbox"/> Signal Violet | <input type="checkbox"/> Sky White | <input type="checkbox"/> Grey Smoke | <input type="checkbox"/> Misty Grey |

FASHION COLOR N/C

- | | | | |
|---------------------------------------|---|---------------------------------------|---|
| <input type="checkbox"/> Candy Blue | <input type="checkbox"/> Candy Purple | <input type="checkbox"/> Candy Red | <input type="checkbox"/> Sparkle Orange |
| <input type="checkbox"/> Glitter Blue | <input type="checkbox"/> Glitter Silver | <input type="checkbox"/> Glitter Gold | <input type="checkbox"/> Black Stardust |
| <input type="checkbox"/> Silver Vein | | | |

SPECTRUM COLOR \$60.00

Additional more than 200 colors available. Please go to www.gunnell-inc.com/PDF/ColorFlyer.pdf to choose.

FRAME WIDTH\$400.00

	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"	
	E2202	<input type="checkbox"/> BT-W2025	25"	\$500.00
	E2202	<input type="checkbox"/> BT-W2026	26"	\$500.00
	E2202	<input type="checkbox"/> BT-W2027	27"	\$500.00
	K0108	<input type="checkbox"/> BT-W2028	28"	\$800.00
	K0108	<input type="checkbox"/> BT-W2029	29"	\$800.00
	K0108	<input type="checkbox"/> BT-W2030	30"	\$1000.00
	K0108	<input type="checkbox"/> BT-W2031	31"	\$1000.00

SEAT DEPTHN/C

	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	
	E2203	<input type="checkbox"/> BT-D2020	20"	\$500.00	
	E2203	<input type="checkbox"/> BT-D2021	21"	\$500.00	
	E2204	<input type="checkbox"/> BT-D2022	22"	\$850.00	
	E2204	<input type="checkbox"/> BT-D2023	23"	\$850.00	
	E2204	<input type="checkbox"/> BT-D2024	24"	\$850.00	
	E2204	<input type="checkbox"/> BT-D2025	25"	\$850.00	
	E2231	<input type="checkbox"/> SEAT PAN		\$100.00	

BACKRESTS OPTIONS

BACK CANE HEIGHTN/C

- (Cane Height is Measured from the Top of the Seat Rail)
- | | | | | | | | | |
|--|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" | <input type="checkbox"/> 21" | <input type="checkbox"/> 22" | <input type="checkbox"/> 23" | <input type="checkbox"/> 24" |
| | <input type="checkbox"/> 25" | <input type="checkbox"/> 26" | <input type="checkbox"/> 27" | | | | | |

K0108 BT-AASH Height & Angle Adjustable Stroller Handle.....\$445.00

√ RECLINE TYPE

- | | | |
|--|--------------------------------|---|
| | <input type="checkbox"/> BT-SR | Static Recline 90° - 130°N/C |
| | | Recline Angle Position |
| | <input type="checkbox"/> 90° | <input type="checkbox"/> 100° <input type="checkbox"/> 110° <input type="checkbox"/> 120° <input type="checkbox"/> 130° |

E1226 BT-75R Dynamic Recline 90° - 180°.....\$945.00

CASTER OPTIONS

- | | | |
|-------|------------------------------------|--|
| | <input type="checkbox"/> BT-515X15 | 5 1/2" X 1 1/2" Composite Casters.....N/C |
| E2214 | <input type="checkbox"/> BT-6X2SPN | 6" X 2" Semi-pneumatic Casters.....\$75.00 |
| E2214 | <input type="checkbox"/> BT-8X2SPN | 8" X 2" Semi-pneumatic Casters.....\$75.00 |
| E1015 | <input type="checkbox"/> BT-FL | Shock Absorbing Casters.....\$900.00 |

DISC BRAKE OR WHEEL LOCK OPTIONS

- | | | |
|-------|---|--|
| | <input type="checkbox"/> BT-WL | Wheel Locks.....N/C |
| E0961 | <input type="checkbox"/> BT-574 | Wheel Lock Extension Handles.....\$42.00/each |
| | <input type="checkbox"/> Disc Brakes (caregiver operated disc brakes) | |
| E2228 | <input type="checkbox"/> BT-DB12 | Disc Brakes with 12" Rear Wheels.....\$975.00 |
| | <input type="checkbox"/> BT-DB20 | Disc Brakes with 20" Rear Wheels\$975.00 |
| | <input type="checkbox"/> BT-DB24 | Disc Brakes with 24" Rear Wheels.....\$975.00 |

WHEEL OPTIONS

AXLES

- | | | |
|-------|---------------------------------|--|
| | <input type="checkbox"/> BT-SA | Standard Axles.....N/C |
| K0108 | <input type="checkbox"/> BT-QRA | Quick Release Axles (only with Wheel Locks)....\$75.00 |

WHEELS

- | | | |
|-------|---------------------------------|--|
| | <input type="checkbox"/> BT-FWD | Reverse Wheel Drive (Rear Wheels in Front).....N/C |
| E2228 | <input type="checkbox"/> BT-12 | 12" Rear Pneumatic Wheels.....N/C |
| E2228 | <input type="checkbox"/> BT-20 | 20" Rear Composite Wheels.....\$325.00 |
| E2228 | <input type="checkbox"/> BT-24 | 24" Rear Composite Wheels.....\$345.00 |
| E2211 | <input type="checkbox"/> BT-PN | Pneumatic Tire.....\$85.00 |
| E2213 | <input type="checkbox"/> BT-FF | Pneumatic Flat Free.....\$65.00 |

WHEEL POSITION N/C

- | | | |
|----------------------------------|---------------------------------|-------------------------------|
| <input type="checkbox"/> Forward | <input type="checkbox"/> Middle | <input type="checkbox"/> Back |
|----------------------------------|---------------------------------|-------------------------------|

SEAT TO FLOOR HEIGHTN/C

- | | | | | | | | |
|--|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| | <input type="checkbox"/> 12" | <input type="checkbox"/> 13" | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" |
| | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" | <input type="checkbox"/> 21" | <input type="checkbox"/> 22" | <input type="checkbox"/> 23" | <input type="checkbox"/> 24" | |

Seat to Floor Height Ranges

		Casters			
		Rear Wheels	5 1/2 x 1 1/2 (BT-515X115)	6 x 2 (BT-6X2PN OR BT-6X2SPN)	8 x 2 (BT-8X2PN OR BT-8X2SPN)
Rear Wheel Drive	12" Wheels (BT-12)	12" - 24"	13" - 24"	15" - 24"	
	20" Wheels (BT-20)	14" - 24"	14" - 24"	15" - 24"	
	24" Wheels (BT-24)	16" - 24"	16" - 24"	16" - 24"	
Reverse Wheel Drive	12" Wheels (BT-12)	12" - 24"	13" - 24"	15" - 24"	
	20" Wheels (BT-20)	14" - 24"	14" - 24"	15" - 24"	
	24" Wheels (BT-24)	16" - 24"	16" - 24"	16" - 24"	

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

F002BRT-1402 Revised May 2014

Rehab Bariatric TNT Adjustable Bariatric Tilt-in-Space Wheelchair

ARMREST OPTIONS (pair)

E0973	<input type="checkbox"/> BT-730	Desk Arms off Uprights.....N/C
E0973	<input type="checkbox"/> BT-705	T-Arms Full Length Height Adjustable.....\$225.00
E0973	<input type="checkbox"/> BT-727	T-Arms Desk Length Height Adjustable.....\$225.00
E0973	<input type="checkbox"/> BT-705DP	Dual Post Full Length Height Adjustable.....\$245.00
E0973	<input type="checkbox"/> BT-727DP	Dual Post Desk Length Height Adjustable.....\$245.00
E0973	<input type="checkbox"/> BT-735	Height & Depth Adjustable Flip up Arms off Uprights <input type="checkbox"/> Full Length.....\$225.00 <input type="checkbox"/> Desk Length.....\$225.00
E0973	<input type="checkbox"/> BT-745	Angle Adjustable Flip up Arms off Uprights, Height and Depth Adjustable <input type="checkbox"/> Full Length.....\$275.00 <input type="checkbox"/> Desk Length.....\$275.00
E0973	<input type="checkbox"/> BT-755	Dynamic Cantilever Armrests.....\$225.00
K0108	<input type="checkbox"/> BT-DPP	Durable Plastic Pads.....\$33.00
K0108	<input type="checkbox"/> BT-441	Gel Pads.....\$120.00

ARMREST HEIGHT OFF SEAT RAIL.....N/C	
<input type="checkbox"/> 8"-12"	<input type="checkbox"/> 8"-16"

FRONT RIGGINGS (pair)

K0052	<input type="checkbox"/> BT-60D	60° Swing Away Leg Rests.....N/C
K0052	<input type="checkbox"/> BT-60	60° Swing Away Lift Off Leg Rests.....\$285.00
K0052	<input type="checkbox"/> BT-70	70° Swing Away Lift Off Leg Rests.....\$285.00
K0052	<input type="checkbox"/> BT-90	90° Swing Away Lift Off Leg Rests.....\$285.00
E0990	<input type="checkbox"/> BT-LSEL	Low-shear Elevating Swing Away Lift Off Leg Rests.....\$549.00
E0990	<input type="checkbox"/> BT-BLC	Bi-Lateral Contracture Footrests.....\$549.00

FOOTPLATES

<input type="checkbox"/> SF	Flip-up Composite Footplates.....N/C	
K0108	<input type="checkbox"/> BT-AASF	Angle Adjustable Single Footplate.....\$165.00
K0040	<input type="checkbox"/> BT-AAF	Angle Adjustable Footplate.....\$195.00/pair
K0108	<input type="checkbox"/> BT-PF	Padded Foot Box.....\$395.00/each
<input type="checkbox"/> BT-CP	Calf Pad.....\$195.00/pair	
<input type="checkbox"/> 440	Ankle Huggers.....\$195.00/pair	

FOOTREST LEGDROP.....N/C	
(Keep Approx. Thickness of Seat Cushion in Mind)	
<input type="checkbox"/> 11"-17"	<input type="checkbox"/> 16"-22"

ADDITIONAL OPTIONS

K0108	<input type="checkbox"/> BT-TP	Transportation Package.....\$250.00
K0108	<input type="checkbox"/> BT-SG	Skirt Guards\$150.00/pair (only with dual post arms options)
K0108	<input type="checkbox"/> BT-LT	Lateral Tilt.....\$3,500.00
K0108	<input type="checkbox"/> BT-TS	Custom Tilt in Space Range.....Quote
<input type="checkbox"/> BT-OAD	One Arm Drive 20" Wheel (choose) <input type="checkbox"/> Left <input type="checkbox"/> Right.....Quote	
<input type="checkbox"/> BT-OAD	One Arm Drive 24" Wheel.... (choose) <input type="checkbox"/> Left <input type="checkbox"/> RightQuote	
K0108	<input type="checkbox"/> BT-AG	Ambulation Guard.....\$130.00 (padded front base frame cross bar)
E0971	<input type="checkbox"/> BT-569T	Anti-tip Device.....\$100.00/pair
<input type="checkbox"/> BT-FSB	Fold Seat to Back Mechanism.....\$250.00	
<input type="checkbox"/> BT-SQR	Seat Quick Removal Mechanism\$250.00	
K0105	<input type="checkbox"/> BT-IVS	IV Pole & Holder.....\$200.00
K0105	<input type="checkbox"/> BT-IVSA	Angle Adjustable IV Pole & Holder.....\$400.00
<input type="checkbox"/> BT-BB	Battery Box.....\$195.00	
E1029	<input type="checkbox"/> BT-VENT	Vent Tray.....\$345.00 (Specify Dimensions) Length____" Width____" Height____"

POWER OPTIONS


E1002	<input type="checkbox"/> BT-PAT	Tilt in Space.....\$4,500.00
K0108	<input type="checkbox"/> BT-PALT	Lateral Tilt.....\$4,500.00
E1003	<input type="checkbox"/> BT-PAR	Recline.....\$4,500.00
E1010	<input type="checkbox"/> BT-PALG	Leg Rest.....\$5,000.00
K0108	<input type="checkbox"/> BT-PALF	Power lift.....Quote

SEATING & POSITIONING


E2607	<input type="checkbox"/> 020S	Solid Seat Base with Standard Hardware.....\$400.00
E2607	<input type="checkbox"/> 020Q	Solid Seat Base with Quick Release Hardware.....\$400.00 Width:____" Depth:____" Padding:____"
E2613	<input type="checkbox"/> 020B	Solid Back Base with Standard Hardware.....\$400.00
E2613	<input type="checkbox"/> 020BQ	Solid Back Base with Quick Release Hardware.....\$400.00 Width:____" Depth:____" Padding:____"

For more options please refer to our "SEATING & POSITIONING" order form at


www.gunnell-inc.com/PDF/SeatingOrderForm.pdf




Approx. Weight 50 lbs/Shipping Weight 75 lbs



Overall Width add 3.5" without handrims and 6.0" with handrims



Overall Height add Seat to Floor Height and Back Cane Height Together



Overall Length Range 22" – 42" Depending on Caster and Wheel Selection
add 3.5" Extra for Anti-tips

FOOTNOTES

- Gunnell will select factory settings if no setting is made
- Dynamic Armrests are recommended with Dynamic Recline Option
- Not all options are compatible with each other. Please contact Gunnell if you have any questions.
- HCPCS codes provided in this form are for reference only. Providers are responsible for determining the appropriate billing codes. Please contact a Medicare professional for further assistance.

**We are a custom company
If you don't see it
Please ask**

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

F002BRT-1402 Revised May 2014