



Rehab RAM

Adjustable Heavy Duty Non-Tilt Wheelchair

Date _____ Account # _____
 Company _____
 Account Contact _____
 Bill to Address _____
 City _____ ST _____ Zip _____
 Tel: _____ Fax _____
 Email: _____

PO # _____
 Ship to _____
 Ship to Attention _____
 Ship to Address _____
 City _____ ST _____ Zip _____
 Tel: _____ Fax _____
 Comments _____

STANDARD FEATURES

- Heavy duty frame and parts (500 lbs Weight Capacity)
- 90° - 130° of adjustable static recline in 10° increments
- 15"-19" seat depth adjustment in 1" increments
- 16"-19" frame width built in adjustment in 1/2" increments
- 15"-19" seat to floor height adjustment in 1" increments
- 17"-27" push handle height adjustment (maintain floor to handle height)
- 8" of forward to wheel base adjustment
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding

K0007 **BASE PRICE**\$2,225.00

FRAME OPTIONS

CLASSIC COLOR N/C

- | | | | |
|--|--|--|--|
| <input type="checkbox"/> Gloss Black | <input type="checkbox"/> Graphite Black | <input type="checkbox"/> Black Wrinkle | <input type="checkbox"/> Lime Green |
| <input type="checkbox"/> Forest Green | <input type="checkbox"/> Cornflower Blue | <input type="checkbox"/> Baby Blue | <input type="checkbox"/> Midnight Blue |
| <input type="checkbox"/> Tangerine | <input type="checkbox"/> Ruby Red | <input type="checkbox"/> Traffic Red | <input type="checkbox"/> Maroon |
| <input type="checkbox"/> Sunny Yellow | <input type="checkbox"/> Bubblegum Pink | <input type="checkbox"/> Raspberry | <input type="checkbox"/> Blue Violet |
| <input type="checkbox"/> Signal Violet | <input type="checkbox"/> Sky White | <input type="checkbox"/> Grey Smoke | <input type="checkbox"/> Misty Grey |

FASHION COLOR N/C

- | | | | |
|---------------------------------------|---|---------------------------------------|---|
| <input type="checkbox"/> Candy Blue | <input type="checkbox"/> Candy Purple | <input type="checkbox"/> Candy Red | <input type="checkbox"/> Sparkle Orange |
| <input type="checkbox"/> Glitter Blue | <input type="checkbox"/> Glitter Silver | <input type="checkbox"/> Glitter Gold | <input type="checkbox"/> Black Stardust |
| <input type="checkbox"/> Silver Vein | | | |

SPECTRUM COLOR \$60.00

Additional more than 200 colors available. Please go to www.gunnell-inc.com/PDF/ColorFlyer.pdf to choose.

FRAME WIDTH N/C

10" 11" 12" 13" 14" 15" 16" 17" 18" 19"

	E2201	RR-W2020	20"\$350.00
	E2201	RR-W2021	21"\$350.00
	E2201	RR-W2022	22"\$350.00
	E2201	RR-W2023	23"\$350.00
	E2201	RR-W2024	24"\$350.00

SEAT DEPTH N/C

11" 12" 13" 14" 15" 16" 17" 18" 19"

	E2203	RR-D2020	20"\$350.00
	E2203	RR-D2021	21"\$350.00
	E2204	RR-D2022	22"\$350.00
	E2204	RR-D2023	23"\$600.00
	E2204	RR-D2024	24"\$600.00
	E2204	RR-D2025	25"\$600.00

E2231 SEAT PAN\$100.00

BACKRESTS OPTIONS

BACK CANE HEIGHT N/C

(Cane Height is measured from the top of the seat rail)
 17" 18" 19" 20" 21" 22" 23" 24"
 25" 26" 27"

K0108 RR-AASH Height & Angle Adjustable Stroller Handl.....\$295.00

RECLINE TYPE

- | | | |
|-------|--------------------------------|--|
| | <input type="checkbox"/> RR-SR | Static Recline 90° - 130°.....N/C |
| | Recline Angle | <input type="checkbox"/> 90° <input type="checkbox"/> 100° <input type="checkbox"/> 110° <input type="checkbox"/> 120° <input type="checkbox"/> 130° |
| E1226 | RR-75R | Dynamic Recline 90° - 180°.....\$695.00 |
| K0108 | RR-75RB | Dynamic Recline with Rocker Back.....\$1295.00 |

CASTER OPTIONS (pair)

- | | | |
|-------|------------------------------------|--|
| | <input type="checkbox"/> RR-515X15 | 5 1/2" X 1 1/2" Composite Casters.....N/C |
| E2219 | RR-3X1 | 3" X 1" Composite Casters.....\$60.00 |
| E2219 | RR-4X1 | 4" X 1" Composite Casters.....\$60.00 |
| E2214 | RR-6X2SPN | 6" X 2" Semi-pneumatic Casters.....\$60.00 |
| E2214 | RR-8X2SPN | 8" X 2" Semi-pneumatic Casters.....\$60.00 |
| E1015 | RR-FL | Shock Absorbing Casters.....\$900.00 |

DISC BRAKE OR WHEEL LOCK OPTIONS (pair)

- | | | |
|-------|--------------------------------------|---|
| | <input type="checkbox"/> RR-WL | Wheel Locks.....N/C |
| E0961 | RR-574 | Wheel Lock Extension Handles.....\$32.00/each |
| | <input type="checkbox"/> DiSc Brakes | (caregiver operated disc brakes) |
| E2228 | RR-DB12 | Disc Brakes with 12" Rear Wheels.....\$975.00 |
| E2228 | RR-DB20 | Disc Brakes with 20" Rear Wheels.....\$975.00 |
| E2228 | RR-DB24 | Disc Brakes with 24" Rear Wheels.....\$975.00 |

WHEEL OPTIONS (pair)

AXLES

- | | | |
|-------|--------------------------------|--|
| | <input type="checkbox"/> RR-SA | Standard Axles.....N/C |
| K0108 | RR-QRA | Quick Release Axles\$75.00
(only work with Wheel Locks) |

WHEELS

- | | | |
|-------|---------------------------------|--|
| | <input type="checkbox"/> RR-FWD | Reverse Wheel DriveN/C
(Rear Wheels in Front) |
| | <input type="checkbox"/> RR-12 | 12" Rear Pneumatic Wheels.....N/C |
| | <input type="checkbox"/> RR-20 | 20" Rear Composite Wheels.....N/C |
| | <input type="checkbox"/> RR-24 | 24" Rear Composite Wheels.....N/C |
| E2211 | RR-PN | Pneumatic Tire.....\$85.00 |
| E2213 | RR-FF | Pneumatic Flat Free.....\$65.00 |

WHEEL POSITION N/C

Forward Middle Back

SEAT TO FLOOR HEIGHT N/C

	<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"
	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"
	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"			

Seat to Floor Height Ranges

	Real Wheels	Casters				
		3x1	4 x 1	5.5 x 1.5	6 x 2	8 x 2
Rear Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"
Reverse Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

F007RAM-1402 Revised May 2014

Rehab RAM

Adjustable Heavy Duty Non-Tilt Wheelchair

ARMREST OPTIONS (pair)

E0973	<input type="checkbox"/> RR-730	Desk Arms off Uprights	N/C
E0973	<input type="checkbox"/> RR-705	T-Arms Full Length Height Adjustable.....	\$225.00
E0973	<input type="checkbox"/> RR-727	T-Arms Desk Length Height Adjustable.....	\$225.00
E0973	<input type="checkbox"/> RR-705DP	Dual Post Full Length Height Adjustable.....	\$245.00
E0973	<input type="checkbox"/> RR-727DP	Dual Post Desk Length Height Adjustable.....	\$245.00
E0973	<input type="checkbox"/> RR-735	Height & Depth Adjustable Flip up Arms off Uprights	
		<input type="checkbox"/> Full Length.....	\$225.00
		<input type="checkbox"/> Desk Length.....	\$225.00
E0973	<input type="checkbox"/> RR-745	Angle Adjustable Flip up Arms off Uprights, Height and Depth Adjustable	
		<input type="checkbox"/> Full Length.....	\$275.00
		<input type="checkbox"/> Desk Length.....	\$275.00
E0973	<input type="checkbox"/> RR-755	Dynamic Cantilever Armrests.....	\$225.00
K0108	<input type="checkbox"/> RR-DPP	Durable Plastic Pads.....	\$33.00
K0108	<input type="checkbox"/> RR-441	Gel Pads.....	\$120.00

ARMREST HEIGHT OFF SEAT RAIL.....		N/C
<input type="checkbox"/> 8"-12"	<input type="checkbox"/> 8"-16"	

FRONT RIGGINGS (pair)

K0052	<input type="checkbox"/> RR-60D	60° Swing Away Leg Rests.....	N/C
K0052	<input type="checkbox"/> RR-60	60° Swing Away Lift Off Leg Rests.....	\$125.00
K0052	<input type="checkbox"/> RR-70	70° Swing Away Lift Off Leg Rests.....	\$125.00
K0052	<input type="checkbox"/> RR-90	90° Swing Away Lift Off Leg Rests.....	\$125.00
E0990	<input type="checkbox"/> RR-LSEL	Low-shear Elevating Swing Away Lift Off Leg Rests.....	\$295.00
E0990	<input type="checkbox"/> RR-BLC	Bi-Lateral Contracture Footrests.....	\$345.00

FOOTPLATES

	<input type="checkbox"/> RR-FP	Flip-up Composite Footplates.....	N/C
K0108	<input type="checkbox"/> RR-AASF	Angle Adjustable Single Footplate.....	\$165.00
K0040	<input type="checkbox"/> RR-AAF	Angle Adjustable Footplate.....	\$145.00/pair
K0108	<input type="checkbox"/> RR-PF	Padded Foot Box.....	\$395.00/each
	<input type="checkbox"/> RR-CP	Calf Pad.....	\$99.00/each
	<input type="checkbox"/> 440	Ankle Huggers.....	\$95.00/each

FOOTREST LEGDROP.....		N/C
(Keep Approx. Thickness of Seat Cushion in Mind)		
<input type="checkbox"/> 11"-17"	<input type="checkbox"/> 16"-22"	

ADDITIONAL OPTIONS

K0108	<input type="checkbox"/> RR-TP	Transportation Package.....	\$250.00
K0108	<input type="checkbox"/> RR-SG	Skirt Guards (with dual post arm option).....	\$150.00/pair
K0108	<input type="checkbox"/> RR-LT	Lateral Tilt.....	\$3,200.00
K0108	<input type="checkbox"/> RR-TS	Tilt in Space.....	\$900.00
K0108	<input type="checkbox"/> RR-TSC	Custom Tilt in Space Range.....	Quote
	<input type="checkbox"/> RR-OAD	One Arm Drive 20" Wheel <input type="checkbox"/> Left <input type="checkbox"/> Right.....	Quote
	<input type="checkbox"/> RR-OAD	One Arm Drive 24" Wheel <input type="checkbox"/> Left <input type="checkbox"/> Right.....	Quote
K0108	<input type="checkbox"/> RR-AG	Ambulation Guard	\$130.00 (padded front base frame cross bar)
E0971	<input type="checkbox"/> RR-569T	Anti-tip Device.....	\$100.00/pair
	<input type="checkbox"/> RRFSB	Fold Seat to Back Mechanism.....	\$250.00
	<input type="checkbox"/> RRSQR	Seat Quick Removal Mechanism	\$250.00
E0078	<input type="checkbox"/> RR652	2" Positioning Seat Belt	\$85.00
	<input type="checkbox"/> RR337	Positioning Belt Pad, Black.....	\$39.00
K0105	<input type="checkbox"/> RRIVS	IV Pole & Holder.....	\$200.00
K0105	<input type="checkbox"/> RRIVSA	Angle Adjustable IV Pole & Holder.....	\$400.00
	<input type="checkbox"/> RRBB	Battery Box.....	\$195.00
E1029	<input type="checkbox"/> RRVENT	Vent Tray.....	\$345.00
		(Specify Dimensions) Length____" Width____" Height____"	

POWER OPTIONS


E1002	<input type="checkbox"/> RR-PAT	Tilt in Space.....	\$4,500.00
K0108	<input type="checkbox"/> RR-PALT	Lateral Tilt.....	\$4,500.00
E1003	<input type="checkbox"/> RR-PAR	Recline.....	\$4,500.00
E1010	<input type="checkbox"/> RR-PALG	Leg Rest.....	\$5,000.00
K0108	<input type="checkbox"/> RR-PALF	Power Lift.....	Quote

SEATING & POSITIONING


E2607	<input type="checkbox"/> 020S	Solid Seat Base with Standard Hardware.....	\$400.00
E2607	<input type="checkbox"/> 020Q	Solid Seat Base with Quick Release Hardware.....	\$400.00
		Width:____" Depth:____" Padding:____"	
E2613	<input type="checkbox"/> 020B	Solid Back Base with Standard Hardware.....	\$400.00
E2613	<input type="checkbox"/> 020BQ	Solid Back Base with Quick Release Hardware.....	\$400.00
		Width:____" Depth:____" Padding:____"	

For more options please refer to our "SEATING & POSITIONING" order form at


www.gunnell-inc.com/PDF/SeatingOrderForm.pdf




Approx. Weight 50 lbs/Shipping Weight 75 lbs



Overall Width add 3.5" without handrims and 6.0" with handrims



Overall Height add Seat to Floor Height and Back Cane Height Together



Overall Length Range 22" - 42" Depending on Caster and Wheel Selection
add 3.5" Extra for Anti-tips

FOOTNOTES

1. Gunnell will select factory settings if no setting is made
2. Dynamic Armrests are recommended with Dynamic Recline Option
3. Not all options are compatible with each other. Please contact Gunnell if you have any questions.
4. HCPCS codes provided in this form are for reference only. Providers are responsible for determining the appropriate billing codes. Please contact a Medicare professional for further assistance.

We are a custom company
If you don't see it
Please ask

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

F007RAM-1402 Revised May 2014