



Rehab TNT Adjustable Tilt-in-Space Wheelchair

Date _____ Account # _____
 Company _____
 Account Contact _____
 Bill to Address _____
 City _____ ST _____ Zip _____
 Tel _____ Fax _____
 Email: _____

PO # _____
 Ship to _____
 Ship to Attention _____
 Ship to Address _____
 City _____ ST _____ Zip _____
 Tel _____ Fax _____
 Comments _____

STANDARD FEATURES

- Heavy duty frame and parts (500 lbs Weight Capacity)
- Infinite adjustment 55° tilt range (default setting -10° to 45°)
- 90° - 130° of adjustable static recline in 10° increments
- 15"-19" seat depth adjustment in 1" increments
- 16"-19" frame width built in adjustment in 1/2" increments
- 15"-19" seat to floor height adjustment in 1" increments
- 17"-27" push handle height adjustment (maintain floor to handle height)
- 8" of forward to wheel base adjustment
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding

E1161 **BASE PRICE**\$2,945.00

FRAME OPTIONS

CLASSIC COLOR N/C

<input type="checkbox"/> Gloss Black	<input type="checkbox"/> Graphite Black	<input type="checkbox"/> Black Wrinkle	<input type="checkbox"/> Lime Green
<input type="checkbox"/> Forest Green	<input type="checkbox"/> Cornflower Blue	<input type="checkbox"/> Baby Blue	<input type="checkbox"/> Midnight Blue
<input type="checkbox"/> Tangerine	<input type="checkbox"/> Ruby Red	<input type="checkbox"/> Traffic Red	<input type="checkbox"/> Maroon
<input type="checkbox"/> Sunny Yellow	<input type="checkbox"/> Bubblegum Pink	<input type="checkbox"/> Raspberry	<input type="checkbox"/> Blue Violet
<input type="checkbox"/> Signal Violet	<input type="checkbox"/> Sky White	<input type="checkbox"/> Grey Smoke	<input type="checkbox"/> Misty Grey

FASHION COLOR N/C

<input type="checkbox"/> Candy Blue	<input type="checkbox"/> Candy Purple	<input type="checkbox"/> Candy Red	<input type="checkbox"/> Sparkle Orange
<input type="checkbox"/> Glitter Blue	<input type="checkbox"/> Glitter Silver	<input type="checkbox"/> Glitter Gold	<input type="checkbox"/> Black Stardust
<input type="checkbox"/> Silver Vein			

SPECTRUM COLOR \$60.00
 Additional more than 200 colors available. Please go to www.gunnell-inc.com/PDF/ColorFlyer.pdf to choose.

FRAME WIDTHN/C

<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
	E2201	RT-W2020	20"	\$350.00				
	E2201	RT-W2021	21"	\$350.00				
	E2201	RT-W2022	22"	\$350.00				
	E2201	RT-W2023	23"	\$350.00				
	E2201	RT-W2024	24"	\$350.00				

SEAT DEPTHN/C

<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
	E2203	RT-D2020	20"	\$350.00			
	E2203	RT-D2021	21"	\$350.00			
	E2204	RT-D2022	22"	\$350.00			
	E2204	RT-D2023	23"	\$600.00			
	E2204	RT-D2024	24"	\$600.00			
	E2204	RT-D2025	25"	\$600.00			
	E2231	SEAT PAN		\$100.00			

BACKRESTS OPTIONS

BACK CANE HEIGHTN/C
 (Cane Height is measured from the top of the seat rail)

<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
<input type="checkbox"/> 25"	<input type="checkbox"/> 26"	<input type="checkbox"/> 27"					

K0108 RT-AASH Height & Angle Adjustable Stroller Handle....\$295.00

RECLINE TYPE

RT-SR Static Recline 90° - 130°.....N/C
 Recline Angle 90° 100° 110° 120° 130°
 E1226 RT-75R Dynamic Recline 90° - 180°.....\$695.00
 K0108 RT-75RB Dynamic Recline with Rocker Back.....\$1295.00



CASTER OPTIONS (pair)

	<input type="checkbox"/> RT-515X15	5 1/2" X 1 1/2" Composite Casters.....	N/C
E2219	<input type="checkbox"/> RT-3X1	3" X 1" Composite Casters.....	\$60.00
E2219	<input type="checkbox"/> RT-4X1	4" X 1" Composite Casters.....	\$60.00
E2214	<input type="checkbox"/> RT-6X2SPN	6" X 2" Semi-pneumatic Casters.....	\$60.00
E2214	<input type="checkbox"/> RT-8X2SPN	8" X 2" Semi-pneumatic Casters.....	\$60.00
E1015	<input type="checkbox"/> RT-FL	Shock Absorbing Casters.....	\$900.00

DISC BRAKE OR WHEEL LOCK OPTIONS (pair)

	<input type="checkbox"/> RT-WL	Wheel Locks.....	N/C
E0961	<input type="checkbox"/> RT-574	Wheel Lock Extension Handles.....	\$32.00/each
		Disc Brakes (caregiver operated disc brakes)	
E2228	<input type="checkbox"/> RT-DB12	Disc Brakes with 12" Rear Wheels.....	\$975.00
E2228	<input type="checkbox"/> RT-DB20	Disc Brakes with 20" Rear Wheels.....	\$975.00
E2228	<input type="checkbox"/> RT-DB24	Disc Brakes with 24" Rear Wheels.....	\$975.00

WHEEL OPTIONS (pair)

AXLES

	<input type="checkbox"/> RT-SA	Standard Axles.....	N/C
K0108	<input type="checkbox"/> RT-QRA	Quick Release Axles.....	\$75.00
		(only work with Wheel Locks)	

WHEELS

	<input type="checkbox"/> RT-FWD	Reverse Wheel Drive.....	N/C
		(Rear Wheels in Front)	
	<input type="checkbox"/> RT-12	12" Rear Pneumatic Wheels.....	N/C
	<input type="checkbox"/> RT-20	20" Rear Composite Wheels.....	N/C
	<input type="checkbox"/> RT-24	24" Rear Composite Wheels.....	N/C
E2211	<input type="checkbox"/> RT-PN	Pneumatic Tire.....	\$85.00
E2213	<input type="checkbox"/> RT-FF	Pneumatic Flat Free.....	\$65.00

WHEEL POSITION N/C

<input type="checkbox"/> Forward	<input type="checkbox"/> Middle	<input type="checkbox"/> Back
----------------------------------	---------------------------------	-------------------------------

SEAT TO FLOOR HEIGHT.....N/C

<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"
<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"
<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"			

Seat to Floor Height Ranges

	Real Wheels	Casters				
		3x1	4 x 1	5.5 x 1.5	6 x 2	8 x 2
Rear Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"
Reverse Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

F008TNT-1402 Revised May 2014

Rehab TNT Adjustable Tilt-in-Space Wheelchair

ARMREST OPTIONS (pair)

E0973	<input type="checkbox"/> RT-730	Desk Arms off Uprights.....N/C
E0973	<input type="checkbox"/> RT-705	T-Arms Full Length Height Adjustable.....\$225.00
E0973	<input type="checkbox"/> RT-727	T-Arms Desk Length Height Adjustable.....\$225.00
E0973	<input type="checkbox"/> RT-705DP	Dual Post Full Length Height Adjustable.....\$245.00
E0973	<input type="checkbox"/> RT-727DP	Dual Post Desk Length Height Adjustable.....\$245.00
E0973	<input type="checkbox"/> RT-735	Height & Depth Adjustable Flip up Arms off Uprights <input type="checkbox"/> Full Length.....\$225.00 <input type="checkbox"/> Desk Length.....\$225.00
E0973	<input type="checkbox"/> RT-745	Angle Adjustable Flip up Arms off Uprights, Height and Depth Adjustable <input type="checkbox"/> Full Length.....\$275.00 <input type="checkbox"/> Desk Length.....\$275.00
E0973	<input type="checkbox"/> RT-755	Dynamic Cantilever Armrests.....\$225.00
K0108	<input type="checkbox"/> RT-DPP	Durable Plastic Pads.....\$33.00
K0108	<input type="checkbox"/> RT-441	Gel Pads.....\$120.00

ARMREST HEIGHT OFF SEAT RAIL.....N/C	
<input type="checkbox"/> 8"-12"	<input type="checkbox"/> 8"-16"

FRONT RIGGINGS (pair)

K0052	<input type="checkbox"/> RT-60D	60° Swing Away Leg Rests.....N/C
K0052	<input type="checkbox"/> RT-60	60° Swing Away Lift Off Leg Rests.....\$125.00
K0052	<input type="checkbox"/> RT-70	70° Swing Away Lift Off Leg Rests.....\$125.00
K0052	<input type="checkbox"/> RT-90	90° Swing Away Lift Off Leg Rests.....\$125.00
E0990	<input type="checkbox"/> RT-LSEL	Low-shear Elevating Swing Away Lift Off Leg Rests.....\$295.00
E0990	<input type="checkbox"/> RT-BLC	Bi-Lateral Contracture Footrests.....\$345.00

FOOTPLATES

<input type="checkbox"/> RT-FP	Flip-up Composite Footplates.....N/C	
K0108	<input type="checkbox"/> RT-AASF	Angle Adjustable Single Footplate.....\$165.00
K0040	<input type="checkbox"/> RT-AAF	Angle Adjustable Footplate.....\$145.00/pair
K0108	<input type="checkbox"/> RT-PF	Padded Foot Box.....\$395.00/each
<input type="checkbox"/> RT-CP	Calf Pad.....\$195.00/pair	
<input type="checkbox"/> 440	Ankle Huggers.....\$195.00/pair	

FOOTREST LEGDROP.....N/C	
(Keep Approx. Thickness of Seat Cushion in Mind)	
<input type="checkbox"/> 11"-17"	<input type="checkbox"/> 16"-22"

ADDITIONAL OPTIONS

K0108	<input type="checkbox"/> RT-TP	Transportation Package.....\$250.00
K0108	<input type="checkbox"/> RT-SG	Skirt Guards\$150.00/pair (only with dual post arms options)
K0108	<input type="checkbox"/> RT-LT	Lateral Tilt.....\$3,200.00
K0108	<input type="checkbox"/> RT-TS	Custom Tilt in Space Range.....Quote
	<input type="checkbox"/> RT-OAD	One Arm Drive 20" Wheel <input type="checkbox"/> Left <input type="checkbox"/> Right.....Quote
	<input type="checkbox"/> RT-OAD	One Arm Drive 24" Wheel <input type="checkbox"/> Left <input type="checkbox"/> Right.....Quote
K0108	<input type="checkbox"/> RT-AG	Ambulation Guard\$130.00 (padded front base frame cross bar)
E0971	<input type="checkbox"/> RT-569T	Anti-tip Device.....\$100.00/pair
	<input type="checkbox"/> RT-FSB	Fold Seat to Back Mechanism.....\$250.00
	<input type="checkbox"/> RT-SQR	Seat Quick Removal Mechanism\$250.00
E0078	<input type="checkbox"/> RT-652	2" Positioning Seat Belt\$85.00
E0978	<input type="checkbox"/> RT-337	Positioning Belt Pad, Black.....\$39.00
K0105	<input type="checkbox"/> RT-IVS	IV Pole & Holder.....\$200.00
K0105	<input type="checkbox"/> RT-IVSA	Angle Adjustable IV Pole & Holder.....\$400.00
	<input type="checkbox"/> RT-BB	Battery Box.....\$195.00
E1029	<input type="checkbox"/> RT-VENT	Vent Tray.....\$345.00 (Specify Dimensions) Length____" Width____" Height____"

POWER OPTIONS


E1002	<input type="checkbox"/> RT-PAT	Tilt in Space.....\$4,500.00
K0108	<input type="checkbox"/> RT-PALT	Lateral Tilt.....\$4,500.00
E1003	<input type="checkbox"/> RT-PAR	Recline.....\$4,500.00
E1010	<input type="checkbox"/> RT-PALG	Leg Rest.....\$5,000.00
K0108	<input type="checkbox"/> RT-PALF	Power LiftQuote

SEATING & POSITIONING


E2607	<input type="checkbox"/> 020S	Solid Seat Base with Standard Hardware.....\$400.00
E2607	<input type="checkbox"/> 020Q	Solid Seat Base with Quick Release Hardware.....\$400.00 Width:____" Depth:____" Padding:____"
E2613	<input type="checkbox"/> 020B	Solid Back Base with Standard Hardware.....\$400.00
E2613	<input type="checkbox"/> 020BQ	Solid Back Base with Quick Release Hardware.....\$400.00 Width:____" Depth:____" Padding:____"

For more options please refer to our "SEATING & POSITIONING" order form at


www.gunnell-inc.com/PDF/SeatingOrderForm.pdf




Approx. Weight 50 lbs/Shipping Weight 75 lbs



Overall Width add 3.5" without handrims and 6.0" with handrims



Overall Height add Seat to Floor Height and Back Cane Height Together



Overall Length Range 22" – 42" Depending on Caster and Wheel Selection

FOOTNOTES

- Gunnell will select factory settings if no setting is made
- Dynamic Armrests are recommended with Dynamic Recline Option
- Not all options are compatible with each other. Please contact Gunnell if you have any questions.
- HCPCS codes provided in this form are for reference only. Providers are responsible for determining the appropriate billing codes. Please contact a Medicare professional for further assistance.

**We are a custom company
If you don't see it
Please ask**

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

F008TNT-1402 Revised May 2014