



Rehab TNT Adjustable Tilt-in-Space Wheelchair

Date _____ Account # _____
 Company _____
 Account Contact _____
 Bill to Address _____
 City _____ ST _____ Zip _____
 Tel _____ Fax _____
 Email: _____

PO # _____
 Ship to _____
 Ship to Attention _____
 Ship to Address _____
 City _____ ST _____ Zip _____
 Tel _____ Fax _____
 Comments _____

STANDARD FEATURES

- Heavy duty frame and parts (500 lbs Weight Capacity)
- Infinite adjustment 55° tilt range (default setting -10° to 45°)
- 90° - 130° of adjustable static recline in 10° increments
- 15"-19" seat depth adjustment in 1" increments
- 16"-19" frame width built in adjustment in 1/2" increments
- 15"-19" seat to floor height adjustment in 1" increments
- 17"-27" push handle height adjustment (maintain floor to handle height)
- 8" of forward to wheel base adjustment
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding

E1161 **BASE PRICE** \$ 3539.00

FRAME OPTIONS

CLASSIC COLOR N/C

- | | | | |
|--|--|--|--|
| <input type="checkbox"/> Gloss Black | <input type="checkbox"/> Graphite Black | <input type="checkbox"/> Black Wrinkle | <input type="checkbox"/> Lime Green |
| <input type="checkbox"/> Forest Green | <input type="checkbox"/> Cornflower Blue | <input type="checkbox"/> Baby Blue | <input type="checkbox"/> Midnight Blue |
| <input type="checkbox"/> Tangerine | <input type="checkbox"/> Ruby Red | <input type="checkbox"/> Traffic Red | <input type="checkbox"/> Maroon |
| <input type="checkbox"/> Sunny Yellow | <input type="checkbox"/> Bubblegum Pink | <input type="checkbox"/> Raspberry | <input type="checkbox"/> Blue Violet |
| <input type="checkbox"/> Signal Violet | <input type="checkbox"/> Sky White | <input type="checkbox"/> Grey Smoke | <input type="checkbox"/> Misty Grey |

FASHION COLOR \$250.00

- | | | | |
|---------------------------------------|---|---------------------------------------|---|
| <input type="checkbox"/> Candy Blue | <input type="checkbox"/> Candy Purple | <input type="checkbox"/> Candy Red | <input type="checkbox"/> Sparkle Orange |
| <input type="checkbox"/> Glitter Blue | <input type="checkbox"/> Glitter Silver | <input type="checkbox"/> Glitter Gold | <input type="checkbox"/> Black Stardust |
| <input type="checkbox"/> Silver Vein | | | |

SPECTRUM COLOR \$60.00

Additional more than 200 colors available. Please go to www.gunnell-inc.com/PDF/ColorFlyer.pdf to choose.

FRAME WIDTH N/C

<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
	E2201	RT-W2020	20"	\$350.00				
	E2201	RT-W2021	21"	\$350.00				
	E2201	RT-W2022	22"	\$350.00				
	E2201	RT-W2023	23"	\$350.00				
	E2201	RT-W2024	24"	\$350.00				

SEAT DEPTH N/C

<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
	E2203	RT-D2020	20"	\$350.00			
	E2203	RT-D2021	21"	\$350.00			
	E2204	RT-D2022	22"	\$350.00			
	E2204	RT-D2023	23"	\$600.00			
	E2204	RT-D2024	24"	\$600.00			
	E2204	RT-D2025	25"	\$600.00			

E2231 SEAT PAN \$100.00

BACKRESTS OPTIONS

BACK CANE HEIGHT N/C

(Cane Height is measured from the top of the seat rail)

- | | | | | | | | |
|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" | <input type="checkbox"/> 21" | <input type="checkbox"/> 22" | <input type="checkbox"/> 23" | <input type="checkbox"/> 24" |
| <input type="checkbox"/> 25" | <input type="checkbox"/> 26" | <input type="checkbox"/> 27" | | | | | |

K0108 RT-AASH Height & Angle Adjustable Stroller Handle....\$595.00

RECLINE TYPE

- | | | |
|--|--|----------|
| <input type="checkbox"/> RT-SR | Static Recline 90° - 130°..... | N/C |
| | Recline Angle <input type="checkbox"/> 90° <input type="checkbox"/> 100° <input type="checkbox"/> 110° <input type="checkbox"/> 120° <input type="checkbox"/> 130° | |
| <input type="checkbox"/> E1226 <input type="checkbox"/> RT-75R | Dynamic Recline 90° - 180°..... | \$695.00 |

CASTER OPTIONS (pair)

- | | | |
|--|--|----------|
| <input type="checkbox"/> RT-515X15 | 5 1/2" X 1 1/2" Composite Casters..... | \$79.00 |
| E2214 <input type="checkbox"/> RT-6X2SPN | 6" X 2" Semi-pneumatic Casters..... | \$89.00 |
| E2214 <input type="checkbox"/> RT-8X2SPN | 8" X 2" Semi-pneumatic Casters..... | \$89.00 |
| <input type="checkbox"/> RT-FL | Shock Absorbing Casters..... | \$900.00 |

DISC BRAKE OR WHEEL LOCK OPTIONS (pair)

- | | | |
|--|---|--------------|
| <input type="checkbox"/> RT-WL | Wheel Locks..... | N/C |
| E0961 <input type="checkbox"/> RT-574 | Wheel Lock Extension Handles..... | \$69.00/each |
| | <input type="checkbox"/> Disc Brakes (caregiver operated disc brakes) | |
| E2228 <input type="checkbox"/> RT-DB12 | Disc Brakes with 12" Rear Wheels..... | \$975.00 |
| E2228 <input type="checkbox"/> RT-DB20 | Disc Brakes with 20" Rear Wheels..... | \$975.00 |
| E2228 <input type="checkbox"/> RT-DB24 | Disc Brakes with 24" Rear Wheels..... | \$975.00 |

WHEEL OPTIONS (pair)

AXLES

- | | | |
|---------------------------------------|------------------------------|---------|
| <input type="checkbox"/> RT-SA | Standard Axles..... | N/C |
| K0108 <input type="checkbox"/> RT-QRA | Quick Release Axles..... | \$89.00 |
| | (only work with Wheel Locks) | |

WHEELS

- | | | |
|--------------------------------------|--|----------------|
| <input type="checkbox"/> RT-FWD | Reverse Wheel Drive ..(Rear wheels in front)..... | N/C |
| | HH-444 Handrim for 20" and 24" Wheels only..... | \$99.00 |
| <input type="checkbox"/> RT-12 | 12" Rear Pneumatic Wheels..... | N/C |
| <input type="checkbox"/> RT-20 | 20" Rear Composite Wheels..... | N/C |
| <input type="checkbox"/> RT-24 | 24" Rear Composite Wheels..... | N/C |
| E2211 <input type="checkbox"/> RT-PN | Pneumatic Tire..... | \$85.00 |
| E2213 <input type="checkbox"/> RT-FF | Pneumatic Flat Free..... | \$89.00 |

WHEEL POSITION

- Forward Middle Back

SEAT TO FLOOR HEIGHT N/C

- | | | | | | |
|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| <input type="checkbox"/> 10" | <input type="checkbox"/> 11" | <input type="checkbox"/> 12" | <input type="checkbox"/> 13" | <input type="checkbox"/> 14" | <input type="checkbox"/> 15" |
| <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" | <input type="checkbox"/> 21" |
| <input type="checkbox"/> 22" | <input type="checkbox"/> 23" | <input type="checkbox"/> 24" | | | |

Seat to Floor Height Ranges

	Real Wheels	Casters				
		5.5 x 1.5	6 x 2	8 x 2		
Rear Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"
Reverse Wheel Drive	12"	10" - 24"	11" - 24"	12" - 24"	13" - 24"	15" - 24"
	20"	14" - 24"	14" - 24"	14" - 24"	14" - 24"	15" - 24"
	24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"	16" - 24"

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

info@gunnell-inc.com Revised Feb 2021

Rehab TNT Adjustable Tilt-in-Space Wheelchair

ARMREST OPTIONS (pair)

E0973	<input type="checkbox"/> RT-730	Desk Arms off Uprights.....	N/C
E0973	<input type="checkbox"/> RT-705	T-Arms Full Length Height Adjustable.....	\$278.00
E0973	<input type="checkbox"/> RT-727	T-Arms Desk Length Height Adjustable.....	\$278.00
E0973	<input type="checkbox"/> RT-705DP	Dual Post Full Length Height Adjustable.....	\$298.00
E0973	<input type="checkbox"/> RT-727DP	Dual Post Desk Length Height Adjustable.....	\$298.00
E0973	<input type="checkbox"/> RT-735	Height & Depth Adjustable Flip up Arms off Uprights	
		<input type="checkbox"/> Full Length.....	\$278.00
		<input type="checkbox"/> Desk Length.....	\$278.00
E0973	<input type="checkbox"/> RT-745	Angle Adjustable Flip up Arms off Uprights, Height and Depth Adjustable	
		<input type="checkbox"/> Full Length.....	\$278.00
		<input type="checkbox"/> Desk Length.....	\$278.00
E0973	<input type="checkbox"/> RT-755	Dynamic Cantilever Armrests.....	\$278.00
K0108	<input type="checkbox"/> RT-DPP	Durable Plastic Pads.....	\$66.00
K0108	<input type="checkbox"/> RT-441	Gel Pads.....	\$369.00

ARMREST HEIGHT OFF SEAT RAIL.....		N/C	
<input type="checkbox"/>	8"-12"	<input type="checkbox"/>	8"-16"

FRONT RIGGINGS (pair)

K0052	<input type="checkbox"/> RT-60D	60° Swing Away Leg Rests.....	\$159.00
K0052	<input type="checkbox"/> RT-60	60° Swing Away Lift Off Leg Rests.....	\$159.00
K0052	<input type="checkbox"/> RT-70	70° Swing Away Lift Off Leg Rests.....	\$159.00
K0052	<input type="checkbox"/> RT-90	90° Swing Away Lift Off Leg Rests.....	\$159.00
E0990	<input type="checkbox"/> RT-LSEL	Low-shear Elevating Swing Away Lift Off Leg Rests.....	\$590.00
E0990	<input type="checkbox"/> RT-BLC	Bi-Lateral Contracture Footrests.....	\$590.00

FOOTPLATES

<input type="checkbox"/>	RT-FP	Flip-up Composite Footplates.....	N/C
K0108	<input type="checkbox"/> RT-AASF	Angle Adjustable Single Footplate.....	\$165.00
K0040	<input type="checkbox"/> RT-AAF	Angle Adjustable Footplate.....	\$165.00/pair
K0108	<input type="checkbox"/> RT-PF	Padded Foot Box.....	\$549.00/each
	<input type="checkbox"/> RT-CP	Calf Pad.....	\$329.00/pair
	<input type="checkbox"/> 440	Ankle Huggers.....	\$218.00/pair

FOOTREST LEGDROP.....		N/C	
(Keep Approx. Thickness of Seat Cushion in Mind)			
<input type="checkbox"/>	11"-17"	<input type="checkbox"/>	16"-22"

ADDITIONAL OPTIONS

K0108	<input type="checkbox"/> RT-TP	Transportation Package.....	\$250.00
K0108	<input type="checkbox"/> RT-SG	Skirt Guards.....	\$295.00/pair (only with dual post arms options)
K0108	<input type="checkbox"/> RT-LT	Lateral Tilt.....	\$3,200.00
K0108	<input type="checkbox"/> RT-TS	Custom Tilt in Space Range.....	Quote
	<input type="checkbox"/> RT-OAD	One Arm Drive 20" Wheel <input type="checkbox"/> Left <input type="checkbox"/> Right.....	Quote
	<input type="checkbox"/> RT-OAD	One Arm Drive 24" Wheel <input type="checkbox"/> Left <input type="checkbox"/> Right.....	Quote
K0108	<input type="checkbox"/> RT-AG	Ambulation Guard.....	\$175.00 (padded front base frame cross bar)
E0971	<input type="checkbox"/> RT-569T	Anti-tip Device.....	\$100.00/pair
	<input type="checkbox"/> RT-FSB	Fold Seat to Back Mechanism.....	\$295.00
	<input type="checkbox"/> RT-SQR	Seat Quick Removal Mechanism.....	\$250.00
E0078	<input type="checkbox"/> RT-652	2" Positioning Seat Belt.....	\$170.00
E0978	<input type="checkbox"/> RT-337	Positioning Belt Pad, Black.....	\$65.00
K0105	<input type="checkbox"/> RT-IVS	IV Pole & Holder.....	\$400.00
K0105	<input type="checkbox"/> RT-IVSA	Angle Adjustable IV Pole & Holder.....	\$475.00
	<input type="checkbox"/> RT-BB	Battery Box.....	\$225.00
E1029	<input type="checkbox"/> RT-VENT	Vent Tray.....	\$345.00
		(Specify Dimensions) Length ____" Width ____" Height ____"	

POWER OPTIONS


E1002	<input type="checkbox"/> RT-PAT	Tilt in Space.....	\$4,500.00
K0108	<input type="checkbox"/> RT-PALT	Lateral Tilt.....	\$4,500.00
E1003	<input type="checkbox"/> RT-PAR	Recline.....	\$4,500.00
E1010	<input type="checkbox"/> RT-PALG	Leg Rest.....	\$5,000.00
K0108	<input type="checkbox"/> RT-PALF	Power Lift.....	Quote

SEATING & POSITIONING


E2607	<input type="checkbox"/> 020S	Solid Seat Base with Standard Hardware.....	\$490.00
E2607	<input type="checkbox"/> 020Q	Solid Seat Base with Quick Release Hardware.....	\$490.00
		Width: ____" Depth: ____" Padding: ____"	
E2613	<input type="checkbox"/> 020B	Solid Back Base with Standard Hardware.....	\$490.00
E2613	<input type="checkbox"/> 020BQ	Solid Back Base with Quick Release Hardware.....	\$490.00
		Width: ____" Depth: ____" Padding: ____"	

For more options please refer to our "SEATING & POSITIONING" order form at


www.gunnell-inc.com/PDF/SeatingOrderForm.pdf




Approx. Weight 50 lbs/Shipping Weight 75 lbs



Overall Width add 3.5" without handrims and 6.0" with handrims



Overall Height add Seat to Floor Height and Back Cane Height Together



Overall Length Range 22" – 42" Depending on Caster and Wheel Selection

FOOTNOTES

- Gunnell will select factory settings if no setting is made
- Dynamic Armrests are recommended with Dynamic Recline Option
- Not all options are compatible with each other. Please contact Gunnell if you have any questions.
- HCPCS codes provided in this form are for reference only. Providers are responsible for determining the appropriate billing codes. Please contact a Medicare professional for further assistance.

**We are a custom company
If you don't see it
Please ask**

Tel: 800-551-0055 Fax: 800-794-5483 www.gunnell-inc.com

Specifications and prices are subject to change without notice.

info@gunnell-inc.com Revised Feb. 2021