

REHAB TNT

Adjustable Tilt-in-Space Wheelchair

Made in the U.S.A.



Date

Account #

Account Contact

Company

First Name

Last Name

Phone Number

Shipping Address

Billing Address

Street Address

Street Address

City

State

City

State

Zip Code

Zip Code

Email

Comments



GUNNELL INC.
MOBILITY SOLUTIONS



PHONE: 800-551-0055



WEBSITE: GUNNELL-INC.COM





Standard Features

Base Price E1161 - \$ 3,679

- Heavy duty frame and parts (500 lbs Weight Capacity)
- Infinite adjustment 55° tilt range (default setting -10° to 45°)
- 90° - 130° of adjustable static recline in 10° increments
- 15"-19" seat depth adjustment in 1" increments
- 16"-19" frame width built in adjustment in 1/2" increments
- 15"-19" seat to floor height adjustment in 1" increments
- 17"-27" push handle height adjustment (maintain floor to handle height)
- 8" of forward to wheel base adjustment
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding

Frame Options

Choose between Standard - Black N/C or Another Color - \$250

Standard - Black N/C



Graphite Black - \$250



Black Wrinkle - \$250



Lime Green - \$250



Forest Green - \$250



Signal Violet - \$250



Baby Blue - \$250



Midnight Blue - \$250



Tangerine - \$250



Ruby Red - \$250



Traffic Red - \$250



Maroon - \$250



Sunny Yellow - \$250



Bubblegum Pink - \$250



Raspberry - \$250



Blue Violet - \$250



Cornflower Blue - \$250



Sky White - \$250




Grey Smoke - \$250



Misty Grey - \$250

Frame Width


16"-19" N/C; 20"-22" \$570; 23"-24" \$800



16"	RT-W2020 20" (E2201) \$570
17"	RT-W2021 21" (E2201) \$570
18"	RT-W2022 22" (E2201) \$570
19"	RT-W2023 23" (E2201) \$800
	RT-W2024 24" (E2201) \$800

Seat Depth

16"-19" N/C; 20"-21" \$620; 23"-25" \$1,030



16"	RT-D2020 20" (E2203) \$620
17"	RT-D2021 21" (E2203) \$620
18"	RT-D2022 22" (E2204) \$1,030
19"	RT-D2023 23" (E2204) \$1,030
	RT-D2024 24" (E2204) \$1,030
	RT-D2025 25" (E2204) \$1,030

Seating & Positioning

Seat Pan (E2231) (pan only, no cushion) \$220

020S Solid Seat Base with Standard Hardware (E2607) \$490

020Q Solid Seat Base with Quick Release Hardware (E2607) \$490

No seat pan, no cushion N/C

Padding, in inches (N/C)

1" 2" 3" 4" 5" 6"

Seating Options

RT-SQR Seat Quick Removal Mechanism (E0971) \$250

RT-652 2" Positioning Seat Belt (E0078) \$70

RT-337 Positioning Belt Pad, Black (E0978) \$170

Back Cane Height (N/C)

Cane Height is measured from the top of the seat rail



17"	21"
18"	22"
19"	23"
20"	24"

Height & Angle Adjustable Stroller

RT-AASH Height & Angle Adjustable Stroller (K0108) \$595

Back Base

020B Solid Back Base with Standard Hardware (E2613) \$490

020BQ Solid Back Base with Quick Release Hardware (E2613) \$490

No Base, No Padding N/C

Base Only, No Padding N/C

Solid - (no back base, no back padding) N/C

Padding Thickness (N/C)

(in inches)

1" 2" 3" 4" 5" 6"

Back Options

RT-IVS IV Pole & Holder (K0105) \$400

RT-IVSA Angle Adjustable IV Pole & Holder (K0105) \$475

Recline Type

Static Recline (90°-130°) or Dynamic Recline (90°-180°)



RT-SR Static Recline N/C

RT-75R Dynamic Recline (E1226) \$1,685

Please Note: Minimum seat depth becomes 17" for Dynamic Recline

Recline Angle (N/C)

90° 100° 110° 120° 130°

Armrest Options (Pair)

NOTE: Dynamic Armrests are recommended with Dynamic Recline Option | Code: (E0973)

RT-730 Desk Arms Off Uprights N/C

RT-705 T-Arms Off Seat Frame Full Length Height Adjustable \$278

RT-727 T-Arms Off Seat Frame Desk Length Height Adjustable \$278

RT-705DP Dual Post Off Seat Frame Full Length Height Adjustable \$298

RT-727DP Dual Post Off Seat Frame Desk Length Height Adjustable \$298

RT-735 Height & Depth Adjustable Flip-Up Arms Off Uprights Full Length \$278

RT-735 Height & Depth Adjustable Flip-Up Arms Off Uprights Desk Length \$278

RT-745 Angle Adjustable Flip-Up Arms Off Uprights, Height & Depth Adjustable Full Length \$278

RT-745 Angle Adjustable Flip-Up Arms Off Uprights, Height & Depth Adjustable Desk Length \$278

RT-755 Dynamic Cantilever Armrests Off Uprights \$278
None

Pad Options

RT-DPP Durable Plastic Pads (K0108) \$66

RT-411 Gel Pads (K0108) \$369

Armrest Height Off Seat Rail (N/C)

8"-12"

8"-16"

N/A

Front Riggings

RT-60D 60° Swing Away Leg Rests (K0052) \$159
RT-60 60° Swing Away Lift Off Leg Rests (K0052) \$159
RT-70 70° Swing Away Lift Off Leg Rests (K0052) \$159
RT-90 90° Swing Away Lift Off Leg Rests (K0052) \$159
RT-LSEL Low-Shear Elevating Swing Away Lift Off Leg Rests (E0990) \$590
RT-BLC Bilateral Contracture Footrests (E0990) \$590
(20" wide chair only)
None

Footplates

RT-FP Flip-Up Composite Footplates N/C
RT-AASF Angle Adjustable Single Footplate (K0108) \$165
RT-AAF Angle Adjustable Footplate (K0040) \$165/pair
RT-PF Padded Foot Box (K0108) \$549/each
None

Footplates Options

RT-CP Calf Pads \$129/pair
FT-240L Ankle Huggers (E0990) \$218/pair
FT-221 Toe Strap - \$49/pair

Footrest Legdrop (N/C)

(keep approx. thickness of seat cushion in mind)

11"-17"
16"-22"
N/A

Caster Options (Pair)

6" x 2" - Default Option | Code: (E2214)

RT-515X15 5½" x 1½" Composite Casters \$79
RT-6X2SPN 6" x 2" Semi-pneumatic Casters \$89
RT-8X2SPN 8" x 2" Semi-pneumatic Casters \$89

Shock Absorbing Casters

RT-FL Shock Absorbing Casters \$900

Wheel Options (Pair)

Axles

RT-SA Standard Axles N/C
RT-QRA Quick Release Axles (K0108) \$89

Please Note: If you selected Disc Brakes, you cannot choose Quick Release Axles.

Wheels

RT-12 12" Rear Pneumatic Wheels
RT-20 20" Rear Composite Wheels
RT-24 24" Rear Composite Wheels

RT-12 12" Rear Pneumatic Wheels Options

RT-FF Pneumatic Flat Free (E2213) \$60
RT-FWD Reverse Wheel Drive (Rear wheels in front) N/C

RT-20 20" Rear Composite Wheels Options

RT-PN Pneumatic Tire (E2211) \$85
RT-FF Pneumatic Flat Free (E2213) \$60
HH-444 Handrim (for 20" & 24" wheels only) \$99

RT-24 24" Rear Composite Wheels Options

RT-PN Pneumatic Tire (E2211) \$85
RT-FF Pneumatic Flat Free (E2213) \$60
HH-444 Handrim (for 20" & 24" wheels only) \$99

One Arm Drive Options

(must select HH-444 Handrim)

RT-OAD One Arm Drive 20" Wheel (Quote)
RT-OAD One Arm Drive 24" Wheel (Quote)

One Arm Drive (N/C)

Left Right

Skirt Guards

RT-SG Skirt Guards (K0108) \$295/pair

Wheel Position (N/C)

Forward Middle Back (Default)

Disc Brake or Wheel Lock Options (Pair)

Disc Brakes

RT-DB12 Disc Brakes with 12" Rear Wheels (E2228) - \$1,460
RT-DB20 Disc Brakes with 20" Rear Wheels (E2228) - \$1,460
RT-DB24 Disc Brakes with 24" Rear Wheels (E2228) - \$1,460

Wheel Locks or Disc Brakes (N/C)

RT-WL Wheel Locks
Disc Brakes (caregiver operated)
None

Wheel Lock Extension Handles

RT-574 Extension for Wheel Locks \$70/each
(Does not apply to disc breaks)

REHAB TNT

Adjustable Tilt-in-Space Wheelchair

Made in the U.S.A.



Seat to Floor Height Ranges:

	Rear Wheels	Front Casters		
		5 x 1.5	6 x 2	8 x 2
Standard Frame	12"	12"-20"	13"-20"	15"-20"
	20"	14"-20"	14"-20"	15"-20"
	24"	16"-20"	16"-20"	15"-20"
Frame with Battery or Vent Tray	12"	16"-20"	16"-20"	16"-20"
	20"	18"-20"	18"-20"	18"-20"
	24"	20"-20"	20"-20"	20"-20"
Standard Frame with Lateral Tilt	12"	16"-24"	17"-24"	19"-24"
	20"	18"-24"	18"-24"	19"-26"
	24"	20"-24"	20"-24"	19"-26"

Seat to Floor Height (N/C)

NOTE: Lateral Tilt or Tilt in Space adds 4" to seat to floor height. Make selection with this value added. Minimum becomes 18"



14"	20"
15"	21"
16"	22"
17"	23"
18"	24"
19"	


Additional Options

RT-TP Transportation Package (K0108) \$360
RT-LB Single Button Actuated Lap Belt (K0108) - \$42
RT-LT Lateral Tilt (K0108) \$3,200
RT-TS Custom Tilt in Space Range (K0108) (Quote)
RT-AG Ambulation Guard (padded front base frame cross bar) (K0108) \$150
RT-569T Anti-tip Device (E0971) \$90/pair
RT-BB Battery Box \$150
RT-VENT Vent Tray (E1029) \$560/ea


Power Options

RT-PAT Tilt in Space (E1002) \$4,500
RT-PALT Lateral Tilt (K0108) \$4,500
RT-PAR Recline (E1003) \$4,500
RT-PALG Leg Rest (E1010) \$5,000
RT-PALF Power Lift (K0108) (Quote)

1. Gunnell will select factory settings if no setting is made.
2. Not all options are compatible with each other. Please contact Gunnell if you have any questions.
3. HCPCS codes provided in this form are for reference only. Providers are responsible for determining the appropriate billing codes. Please contact a Medicare professional for further assistance.
4. 4. Approx. Weight 50 lbs/Shipping Weight 75 lbs



Overall Width add 3.5" without handrims and 6.0" with handrims



Overall Height add Seat to Floor Height and Back Cane Height Together



Overall Length Range 22" – 42" Depending on Caster and Wheel Selection



GUNNELL INC.
MOBILITY SOLUTIONS



PHONE: 800-551-0055



WEBSITE: GUNNELL-INC.COM

