

REHAB KIDSTER

Adjustable Pediatric Tilt-in-Space Wheelchair

Made in the U.S.A.



Date

Account #

Account Contact

Company

First Name

Last Name

Phone Number

Shipping Address

Billing Address

Street Address

Street Address

City

State

City

State

Zip Code

Zip Code

Email

Comments



GUNNELL INC.
MOBILITY SOLUTIONS



PHONE: 800-551-0055



WEBSITE: GUNNELL-INC.COM





Standard Features

Base Price E1161 - \$ 2,989

- Heavy duty frame and parts (500 lbs Weight Capacity)
- Infinite adjustment 55° tilt range (default setting -10° to 45°)
- 90° - 130° of adjustable static recline in 10° increments
- 15"-19" seat depth adjustment in 1" increments
- 16"-19" frame width built in adjustment in 1/2" increments
- 15"-19" seat to floor height adjustment in 1" increments
- 17"-24" push handle height adjustment (maintain floor to handle height)
- 8" of forward to wheel base adjustment
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding

Frame Options

Choose between Standard - Black N/C or Another Color - \$250

Standard - Black N/C



Graphite Black - \$250



Black Wrinkle - \$250



Lime Green - \$250



Forest Green - \$250



Signal Violet - \$250



Baby Blue - \$250



Midnight Blue - \$250



Tangerine - \$250



Ruby Red - \$250



Traffic Red - \$250



Maroon - \$250



Sunny Yellow - \$250



Bubblegum Pink - \$250



Raspberry - \$250



Blue Violet - \$250



Cornflower Blue - \$250



Sky White - \$250




Grey Smoke - \$250



Misty Grey - \$250

Frame Width


10"-13" - Custom; 14"-17" - N/C



10" - Custom	14"
11" - Custom	15"
12" - Custom	16"
13" - Custom	17"

Seat Depth

11"-12" - Custom; 13"-17" - N/C



11" - Custom	13"
12" - Custom	14"
	15"
	16"
	17"

Seating & Positioning

Seat Pan (E2231) (pan only, no cushion) - \$349

020S Solid Seat Base with Standard Hardware (E2607) - \$568

020Q Solid Seat Base with Quick Release Hardware (E2607) - \$568

No seat pan, no cushion - N/C

Padding, in inches (N/C)

1" 2" 3" 4" 5" 6"

Seating Options

K-SQR Seat Quick Removal Mechanism (E0971) - \$250

K-652 2" Positioning Seat Belt (E0078) - \$84

K-337 Positioning Belt Pad, Black (E0978) - \$65

Back Cane Height

Cane Height is measured from the top of the seat rail



17"	21"	25" - Custom
18"	22"	26" - Custom
19"	23"	27" - Custom
20"	24"	

17"-24" - N/C; 25"-27" - Custom

Height & Angle Adjustable Stroller

K-AASH Height & Angle Adjustable Stroller (K0108) - \$595

Back Base

020B Solid Back Base with Standard Hardware (E2613) - \$770

020BQ Solid Back Base with Quick Release Hardware (E2613) - \$770

No Base, No Padding - N/C

Base Only, No Padding - N/C

Solid - (no back base, no back padding) - N/C

Padding Thickness (N/C)

(in inches)

1" 2" 3" 4" 5" 6"

Back Options

K-IVS IV Pole & Holder (K0105) - \$226

K-IVSA Angle Adjustable IV Pole & Holder (K0105) - \$226

Recline Type

Static Recline (90°-130°) or Dynamic Recline (90°-180°)



K-SR Static Recline - N/C

K-75R Dynamic Recline (E1226) - \$1,276

Please Note: Minimum seat depth becomes 17" for Dynamic Recline

Recline Angle (N/C)

90° 100° 110° 120° 130°

Armrest Options (Pair)

NOTE: Dynamic Armrests are recommended with Dynamic Recline Option | Code: (E0973)

K-730 Desk Arms Off Uprights - N/C

K-735 Height & Depth Adjustable Flip-Up Arms Off Uprights Full Length - \$182

K-735 Height & Depth Adjustable Flip-Up Arms Off Uprights Desk Length - \$182

K-745 Angle Adjustable Flip-Up Arms Off Uprights, Height & Depth Adjustable Full Length - \$182

K-745 Angle Adjustable Flip-Up Arms Off Uprights, Height & Depth Adjustable Desk Length - \$182

K-755 Dynamic Cantilever Armrests Off Uprights - \$182

None

Pad Options

K-DPP Durable Plastic Pads (K0108) - \$182

K-411 Gel Pads (K0108) - \$182

Armrest Height Off Seat Rail (N/C)

6"-14"

Front Riggings

K-60D 60° Swing Away Leg Rests (K0052) - \$182
K-60 60° Swing Away Lift Off Leg Rests (K0052) - \$182
K-70 70° Swing Away Lift Off Leg Rests (K0052) - \$182
K-90 90° Swing Away Lift Off Leg Rests (K0052) - \$182
K-LSEL Low-Shear Elevating Swing Away Lift Off Leg Rests (E0990) - \$259
K-BLC Bilateral Contracture Footrests (E0990) - \$259
(20" wide chair only)
None

Footplates

K-FP Flip-Up Composite Footplates - N/C
K-AASF Angle Adjustable Single Footplate (K0108) - \$165
K-AAF Angle Adjustable Footplate (K0040) - \$146/pair
K-PF Padded Foot Box (K0108) - \$549/each
None - N/C

Footplates Options

K-CP Calf Pads - \$329/pair
FT-240L Ankle Huggers (E0990) - \$218/pair
FT-221 Toe Strap - \$49/pair

Footrest Legdrop (N/C)

(keep approx. thickness of seat cushion in mind)

3"-11"
12"-17"
N/A

Caster Options (Pair)

6" x 2" - Default Option | Code: (E2214)

K-515X15 5½" x 1½" Composite Casters - \$79
K-6X2SPN 6" x 2" Semi-pneumatic Casters - \$89
K-8X2SPN 8" x 2" Semi-pneumatic Casters - \$89

Shock Absorbing Casters

K-FL Shock Absorbing Casters - \$260

Wheel Options (Pair)

Axles

K-SA Standard Axles - N/C
K-QRA Quick Release Axles (K0108) - \$89

Please Note: If you selected Disc Brakes, you cannot choose Quick Release Axles.

Wheels

K-12 12" Rear Pneumatic Wheels
K-20 20" Rear Composite Wheels
K-24 24" Rear Composite Wheels

K-12 12" Rear Pneumatic Wheels Options

K-FF Pneumatic Flat Free (E2213) - \$69
RT-FWD Reverse Wheel Drive (Rear wheels in front) - N/C

K-20 20" Rear Composite Wheels Options

RT-PN Pneumatic Tire (E2211) - \$91
RT-FF Pneumatic Flat Free (E2213) - \$69
HH-444 Handrim (for 20" & 24" wheels only) - \$99

K-24 24" Rear Composite Wheels Options

RT-PN Pneumatic Tire (E2211) - \$91
RT-FF Pneumatic Flat Free (E2213) - \$69
HH-444 Handrim (for 20" & 24" wheels only) - \$99

One Arm Drive Options

(must select HH-444 Handrim)

K-OAD One Arm Drive 20" Wheel - Quote
K-OAD One Arm Drive 24" Wheel - Quote

One Arm Drive (N/C)

Left Right

Skirt Guards

RT-SG Skirt Guards (K0108) - \$295/pair

Wheel Position (N/C)

Forward Middle Back (Default)

Disc Brake or Wheel Lock Options (Pair)

Disc Brakes - \$2,128

K-DB12 Disc Brakes with 12" Rear Wheels (E2228)
K-DB20 Disc Brakes with 20" Rear Wheels (E2228)
K-DB24 Disc Brakes with 24" Rear Wheels (E2228)

Wheel Locks or Disc Brakes (N/C)

K-WL Wheel Locks
Disc Brakes (caregiver operated)
None

Wheel Lock Extension Handles

K-574 Extension for Wheel Locks - \$70/each
(Does not apply to disc breaks)

REHAB KIDSTER

Adjustable Pediatric Tilt-in-Space Wheelchair

Made in the U.S.A.



Seat to Floor Height Ranges:

	Rear Wheels	Front Casters		
		5 x 1.5	6 x 2	8 x 2
Standard Frame	12"	12"-20"	13"-20"	15"-20"
	20"	14"-20"	14"-20"	15"-20"
	24"	16"-20"	16"-20"	15"-20"
Frame with Battery or Vent Tray	12"	16"-20"	16"-20"	16"-20"
	20"	18"-20"	18"-20"	18"-20"
	24"	20"-20"	20"-20"	20"-20"
Standard Frame with Lateral Tilt	12"	16"-24"	17"-24"	19"-24"
	20"	18"-24"	18"-24"	19"-26"
	24"	20"-24"	20"-24"	19"-26"

Seat to Floor Height

NOTE: Lateral Tilt or Tilt in Space adds 4" to seat to floor height. Make selection with this value added. Minimum becomes 18"



10" - Custom 14"
11" - Custom 15"
12" 16"
13"

10"-11" - Custom; 12"-16" - N/C


Additional Options

K-TP Transportation Package (K0108) - \$369
K-LB Single Button Actuated Lap Belt (K0108) - \$42
K-LT Lateral Tilt (K0108) - \$3,200
K-TS Custom Tilt in Space Range (K0108) - Quote
K-AG Ambulation Guard (padded front base frame cross bar) (K0108) - \$175
K-569T Anti-tip Device (E0971) - \$98/pair
K-BB Battery Box - \$226
K-VENT Vent Tray (E1029) - \$724/ea


Power Options

K-PAT Tilt in Space (E1002) - \$4,500
K-PALT Lateral Tilt (K0108) - \$4,500
K-PAR Recline (E1003) - \$4,500
K-PALG Leg Rest (E1010) - \$5,000
K-PALF Power Lift (K0108) - Quote

1. Gunnell will select factory settings if no setting is made.
2. Not all options are compatible with each other. Please contact Gunnell if you have any questions.
3. HCPCS codes provided in this form are for reference only. Providers are responsible for determining the appropriate billing codes. Please contact a Medicare professional for further assistance.
4. Approx. Weight 50 lbs/Shipping Weight 75 lbs



Overall Width add 3.5" without handrims and 6.0" with handrims



Overall Height add Seat to Floor Height and Back Cane Height Together



Overall Length Range 22" - 42" Depending on Caster and Wheel Selection



GUNNELL INC.
MOBILITY SOLUTIONS



PHONE: 800-551-0055



WEBSITE: GUNNELL-INC.COM

