

Date

Account #

Account Contact

Company

First Name

Last Name

Phone Number

Shipping Address

Billing Address

Street Address

Street Address

City

State

City

State

Zip Code

Zip Code

Email

Comments

Note: To avoid customer confusion, we are separating the head and foot sides from the left and right sides. Head and foot barriers are now referred to as "Rails" and are selected separately from the sides. "Sides" now only refers to the left and/or right barriers, if selected.





Standard Features

Base Price E1161 - \$ 5,379.75

- Heavy duty frame and parts (400 lbs Weight Capacity)
- Infinite adjustment 45° tilt range (22.5° anterior/22.5° posterior)
- 8" Fixed Head and Foot Rails
- 8" Casters
- 20" Rear Wheels
- Height Adjustable Stroller Handle
- 5 year warranty on all Gunnell designed products
- 90 day warranty on all padding
- Anti-tips
- Manual Brakes

Frame Options

Choose between Standard - Black N/C or Another Color - \$250

Standard - Black N/C



Graphite Black - \$250



Black Wrinkle - \$250



Lime Green - \$250



Forest Green - \$250



Signal Violet - \$250



Baby Blue - \$250



Midnight Blue - \$250



Tangerine - \$250



Ruby Red - \$250



Traffic Red - \$250



Maroon - \$250



Sunny Yellow - \$250



Bubblegum Pink - \$250



Raspberry - \$250



Blue Violet - \$250



Cornflower Blue - \$250



Sky White - \$250



Grey Smoke - \$250



Misty Grey - \$250



Platform Width

(Dimension IDW)

PW 27"-Standard Platform Width - \$1,113.75
 PPCM-001-Custom Platform Width - \$1,410.75

Platform Length

(Specify Dimensions) (Dimensions IDL)

PPCM-002 - Platform Length in inches

Back (IDLB) Seat (IDLS) Leg (IDLL)

Overall Total (IDLB+IDLS+IDLL)

PPCM-002 If the Overall Total is 60" or under - \$1,113.75
 PPCM-002A If the Overall Total is greater than 60" - 70" - \$1,410.75
 PPCM-002B If the Overall Total is greater than 71" - 80" - \$1,559.25

Platform to Floor Height

(Dimension H)

PFH 30" Platform to Floor Height - \$1,113.75
 PPH-001 Custom Platform to Floor Height - \$1,687.50

Hip to Back Angle

HBA-006 Static Hip to Back Angle 180° to 135° - \$324
 HBA-007 Static Hip to Back Angle 180° to 90° - \$384.75
 PPO-008- Static Hip to Back Angle 180° to 135° - \$890
 PPO-010- Static Hip to Back Angle 180° to 90° - \$1,350

Knee Angle

PPKA-006 Static Knee Angle 180° to 135° - \$175.50
 PPKA-007 Static Knee Angle 180° to 90° - \$236.25
 PPO-008- Dynamic Knee Angle 180° to 135° - \$1,485
 PPO-010- Dynamic Knee Angle 180° to 90° - \$1,755

Platform Rails

PPO-001 8" Fixed Head Rail - N/C
 PPO-001 8" Fixed Foot Rail - N/C
 PPCM-005 Custom Height Rails
 (MUST SPECIFY) Desired Height # of sides x - \$75



Platform Sides

(Dimension SH)

PPO-001F 8" Fixed Sides
 # of sides x - \$337.50
 PPO-002 Flip Down Hinged Sides with Dual Locking Clamps
 # of sides x - \$506.25
 PPO-002A 8" Angle Adjustable Positive Locking Sides
 # of sides x - \$837
 PPO-002R 8" Removable Platform Sides
 # of sides x - \$648
 PPCM-002RA 8" Removable Angle Adjustable Positive Locking Sides
 # of sides x - \$1,053
 PPCM-005 Custom Height Sides
 (MUST SPECIFY) Desired Height # of sides x - \$75

Platform Tilt in Space

PTIS Anterior (22.5°)/Posterior(22.5°) Tilt in Space - N/C
 PPCM-006 Custom Anterior(A)/Posterior(P) Tilt in Space - \$1,080

45°A / 0°P	40°A / 5°P	35°A / 10°P	30°A / 15°P
25°A / 20°P	20°A / 25°P	15°A / 30°P	
10°A / 35°P	5°A / 40°P	0°A / 45°P	

Wheel Options

PPRW 20" Rear Wheels - N/C
 PPO-004 24" Wheels (only available on heights 26"+) - \$472.50
 PPO-005 Flat Free Wheel Option - \$168.75

Stroller Handle Options

PSH Height Adjustable Stroller Handle - N/C
 PPO-003 Height & Angle Adjustable Stroller Handle - \$808.65

Custom Options

PPCO-003 Lateral Tilt - \$5,130
 PPCO-TP Transportation Package - \$1,687.50
 PPCO-005 Tapered Corners - \$506.25
 PPCO-012 Roller Bumpers - \$506.25
 PPCO-008 IV-Pole - \$540
 PPCO-001 Vent Tray - \$538.65

(Specify Dimensions in inches)

Length Width Heights

Power Options

PPCO-PATS Tilt in Space - \$6,075
 PPCO-PALT Lateral Tilt - \$6,075
 PPCO-PAR Recline - \$6,075
 PPCO-PALG Legrest - \$6,750
 PPCO-PALF Power Lift - \$12,757.50

Back Platform Base Padding Options

PPBP-001 1" Poly Base - \$290.25
PPBP-002 1" Sunmate Base Pad (Specify Density) - \$425.25
Soft Medium
PPBP-003 2" Poly Base - \$317.25
PPBP-004 2" Sunmate Base (Specify Density) - \$587.25
Soft Medium
PPBP-005 Vinyl Padding Cover with Zipper - \$255.15
PPBP-006 Dartex Padding Cover with Zipper - \$472.50
PPBP-007 Rubatex Padding Cover with Zipper - \$592.65
PPBP-009 Custom Base Padding - Quote

Seat Platform Base Padding Options

PPSP-001 1" Poly Base - \$290.25
PPSP-002 1" Sunmate Base Pad (Specify Density) - \$425.25
Soft Medium
PPSP-003 2" Poly Base - \$303.75
PPSP-004 2" Sunmate Base (Specify Density) - \$587.25
Soft Medium
PPSP-005 Vinyl Padding Cover with Zipper - \$255.15
PPSP-006 Dartex Padding Cover with Zipper - \$472.50
PPSP-007 Rubatex Padding Cover with Zipper - \$592.65
PPSP-009 Custom Base Padding - Quote

Leg Platform Base Padding Options

PPFP-001 1" Poly Base - \$290.25
PPFP-002 1" Sunmate Base Pad (Specify Density) - \$425.25
Soft Medium
PPFP-003 2" Poly Base - \$303.75
PPFP-004 2" Sunmate Base (Specify Density) - \$587.25
Soft Medium
PPFP-005 Vinyl Padding Cover with Zipper - \$255.15
PPFP-006 Dartex Padding Cover with Zipper - \$472.50
PPFP-007 Rubatex Padding Cover with Zipper - \$592.65
PPFP-009 Custom Base Padding - Quote

Rail Padding Options

PPSO-001 1/2" Poly Rail Pads (each)
of sides x - \$195.75
PPSO-002 1/2" Sunmate Rail Pads (Specify Density & Quantity)
Soft Medium
of sides x - \$303.75
PPSO-003 1" Poly Rail Pads (each)
of sides x - \$263.25
PPSO-004 1" Sunmate Rail Pads (Specify Density & Quantity)
Soft Medium
of sides x - \$438.75
PPSO-008 Custom Base Padding - Quote

Rail Cover Options

PPSO-005 Vinyl Padding Cover with Zipper (each) - \$162
PPSO-006 Dartex Padding Cover with Zipper - \$324 # of sides

Side Padding Options

PPSO-001 1/2" Poly Side Pads (each)
of sides x - \$195.75
PPSO-002 1/2" Sunmate Side Pads (Specify Density & Quantity)
Soft Medium
of sides x - \$303.75
PPSO-003 1" Poly Side Pads (each)
of sides x - \$263.25
PPSO-004 1" Sunmate Side Pads (Specify Density & Quantity)
Soft Medium
of sides x - \$438.75
PPSO-008 Custom Base Padding - Quote

Side Cover Options

PPSO-005 Vinyl Padding Cover with Zipper (each) - \$162
PPSO-006 Dartex Padding Cover with Zipper - \$324 # of sides

Footnotes

1. Gunnell will select factory settings if no setting is made.
2. Inside Width (IDW: B, S, or L): The measurement between side pads, side ABS, or width of top platform, whichever is smallest.
3. Inside Length (IDL: B, S, or L): The measurement between side pads, side ABS, or width of top platform, whichever is smallest.
4. Outside Width (ODW: B, S, or L): Inside Width, IDW, plus the following measurement whenever applicable:
(i) 2" for fixed sides;
(ii) 3" for hinged flip down sides;
(iii) 4" for removable sides when sides not removed;
(iv) 5" for angle adjustable locking sides; and
(v) padding thickness
5. Outside Width is the Inside Width for removable sides when sides removed.
6. Platform to Floor Height: measurement from the ground to the top of the platform below the ABS.
7. Maximum 8 sides/rails, left side 1, left side 2, left side 3, right side 4, right side 5, right side 6, head rail, and/or foot rail.
8. Side/Rail Height: measurement from the top of sides to the top of steel platform, ABS, or base pads, whichever is the smallest.
9. Factory setting of 8" Side/Rail Height is measured from top of steel platform to the top of sides/rails.
10. Seat and Leg can be combined as one platform without hinge point in between. Cross out leg platform, initial next to it, write full length in seat position, and cross out knee angle and initial next to it.
11. Static angle adjustments require tools
12. Dynamic angle adjustments are released by a cable assembly and require no tools.

REHAB MAC

Positioning Platform Option

Made in the U.S.A.



Finished Product Dimensions

Length of Back (inside pad to hinge) (*in inches*) IDLB:

Length of Seat (hip angle hinge to knee angle hinge) (*in inches*) IDLS:

Length of Leg Rest (hinge to inside pad) (*in inches*) IDLL:

Width of Back (inside of pads) (*in inches*) IDWB:


Width of Seat (inside of pads) (*in inches*) IDWS:

Width of Leg Rest (inside of pads) (*in inches*) IDWL:

Height of Side (see notes 7,8,9 above) (*in inches*) SH:

Platform to Floor Height (see note 6 above) (*in inches*) H:

Maximum allowable cart width (maximum allowable cart width for going through doors) (*in inches*) ODW:



Platform Sides/Rails

Number of Sides (including head and foot):

- | | |
|------------------|-------------------|
| Head Rail | Foot Rail |
| Back Left Side 1 | Back Right Side 4 |
| Seat Left Side 2 | Seat Right Side 5 |
| Leg Left Side 3 | Leg Right Side 6 |

Options

- Base pad thickness (*in inches*)
- Side pad thickness (*in inches*)



GUNNELL INC.
MOBILITY SOLUTIONS



PHONE: 800-551-0055



WEBSITE: GUNNELL-INC.COM

