

SELF PROPELLING GURNEY OPTION

Rehab MAC

Made in the U.S.A.



Date

Account #

Account Contact

Company

First Name

Last Name

Phone Number

Shipping Address

Billing Address

Street Address

Street Address

City

State

City

State

Zip Code

Zip Code

Email

Comments

Note: To avoid customer confusion, we are separating the head and foot sides from the left and right sides. Head and foot barriers are now referred to as "Rails" and are selected separately from the sides. "Sides" now only refers to the left and/or right barriers, if selected.



GUNNELL INC.
MOBILITY SOLUTIONS



PHONE: 800-551-0055



WEBSITE: GUNNELL-INC.COM





Standard Features

Base Price SPG-001 - \$4,353.75

- Heavy duty frame and parts (400 lbs Weight Capacity)
- 8" Casters
- 20" Rear Wheels
- Height Adjustable Stroller Handle
- 5 year warranty on all Gunnell designed products.
- 90 Days on All Padding
- Anti-tips
- Manual Brakes

Frame Options

Choose between Standard - Black or Another Color

Standard - Black



Graphite Black - \$250



Black Wrinkle - \$250



Lime Green - \$250



Forest Green - \$250



Signal Violet - \$250



Baby Blue - \$250



Midnight Blue - \$250



Tangerine - \$250



Ruby Red - \$250



Traffic Red - \$250



Maroon - \$250



Sunny Yellow - \$250



Bubblegum Pink - \$250



Raspberry - \$250



Blue Violet - \$250



Cornflower Blue - \$250



Sky White - \$250



Grey Smoke - \$250



Misty Grey - \$250



Platform Width

(Dimension IDW)

SPW 22"-Standard Platform Width - \$1,012.50

SPCW-001-Custom Platform Width - \$1,282.50 *In Inches*

Platform Length

(Specify Dimensions) (Dimensions IDL)

SPL 70"-Standard Platform Length - \$1,012.50

SPCL-001-Custom Platform Width - \$1,282.50 *In Inches*

Platform to Floor Height

(Dimension H)

SPFH20 18" Platform to Floor Height (with 20" wheels) - \$1,012.50

SPFH24 22" Platform to Floor Height (with 24" wheels) - \$1,147.50

SPFH-001 Custom Platform to Floor Height - \$1,552.50 *In Inches*

Platform Rails

SPS-001 Foot Rail - N/C

SPS-001F 8" Fixed Head Rail - \$668

SPS-002 8" Flip Foot Rail - \$1,046

SPS-002 8" Flip Head Rail - \$1,046

Platform Sides - Optional

(Dimension SH)

SPS-001F 8" Fixed Sides - \$668.25/side

of sides

SPS-002 8" Flip Down Hinged Sides - \$1,046.25/side

of sides

SPS-003 8" Adjustable Angle Positive Locking Sides - \$1,485/side

of sides

SPS-002R 8" Removable Platform Sides - \$1,188/side

of sides

SPS-002RA 8" Removable Angle Adjustable Positive Locking Sides - \$1,728/side

of sides

SPS -005 Custom Height Sides - \$168.75/side

(MUST SPECIFY) Desired Height # of sides

Platform Base Padding Options

SPBP-001 1" Poly Base - \$472.50

SPBP-002 1" Sunmate Base Pad (Specify Density) - \$810

Soft Medium

SPBP-003 2" Poly Base - \$607.50

SPBP-004 2" Sunmate Base (Specify Density) - \$1,215

Soft Medium

SPBP-009 Custom Base Padding - Quote

Platform Cover Options

SPBP-005 Vinyl Padding Cover with Zipper - \$519.75

SPBP-006 Dartex Padding Cover with Zipper - \$945

SPBP-007 Rubatex Padding Cover with Zipper - \$1,296

Rail Padding Options - Optional

SPSP-001 1/2" Poly Rail Pads - \$391.50 each *QTY*

SPSP-002 1/2" Sunmate Rail Pads (specify density) - \$607.50 each

Soft Medium *QTY*

SPSP-003 1" Poly Rail Pads - \$526.50 each *QTY*

SPSP-004 1" Sunmate Rail Pads (specify density) - \$877.50 each

Soft Medium *QTY*

Side Padding Options - Optional

SPSP-001 1/2" Poly Side Pads - \$391.50 each *QTY*

SPSP-002 1/2" Sunmate Side Pads (specify density) - \$607.50 each

Soft Medium *QTY*

SPSP-003 1" Poly Side Pads - \$526.50 each *QTY*

SPSP-004 1" Sunmate Side Pads (specify density) - \$877.50 each

Soft Medium *QTY*

Rail Cover Options - Optional

SPSP-005 Vinyl Padding Cover with Zipper - \$162 each *QTY*

SPSP-006 Darlex Padding Cover with Zipper - \$648 each *QTY*

SPSP-007 Rubatex Padding Cover with Zipper \$592.65 each *QTY*

SPSP-008 Custom Rail Padding - Quote

Side Cover Options - Optional

SPSP-005 Vinyl Padding Cover with Zipper - \$162 each *QTY*

SPSP-006 Darlex Padding Cover with Zipper - \$648 each *QTY*

SPSP-007 Rubatex Padding Cover with Zipper \$592.65 each *QTY*

SPSP-008 Custom Side Padding - Quote

Wheel Options

SPW20 20" Rear Wheels - N/C

SPR24 24" Rear Wheels - \$472.50

Optional Add-On: SP-FF Flat Free Wheel Option - \$168.75

Stroller Handle Options

SPSH Height Adjustable Stroller Handle - N/C

SPO-003 Height & Angle Adjustable Stroller Handle - \$808.65

Custom Options

Select any of the following:

SPCO-TS Tilt in Space - \$5,130

SPCO-003 Lateral Tilt - \$5,130

SPCO-TP Transportation Package - \$1,687.50

SPCO-001 Tapered Corners - \$506.25

SPCO-012 Roller Bumpers - \$506.25

SPCO-002 IV-Pole - \$540

SPCO-013 Anti Tip - \$100

SPCO-019 Vent Tray - \$538.65

Power Options

Select any of the following:

SPCO-PATS Tilt in Space - \$6,075


SPCO-PALT Lateral Tilt - \$6,075

SPCO-PALF Power Lift - \$12,757.50

SELF PROPELLING GURNEY OPTION

Rehab MAC

Made in the U.S.A.



Finished Product Dimensions

Inside Length IDL:

Inside Width IDW:

Outside Length ODL:

Outside Width ODW:

Height of Side (see notes 7,8,9 above) SH:

Platform to Floor Height (see note 6 above) H:

Maximum allowable cart width (maximum allowable cart width for going through doors) ODW:



Platform Sides/Rails

Number of Sides (including head and foot):

Head Rail

Foot Rail

Left Side

Right Side

Options

Base pad thickness (*in inches*)

Side pad thickness (*in inches*)

Footnotes

1. Gunnell will select factory settings if no setting is made

2. Inside Width (IDW): The measurement between side pads, side ABS, or width of top platform, whichever is smallest

3. Inside Length (IDL): The measurement between side pads, side ABS, or width of top platform, whichever is smallest

4. Outside Width (ODW): Inside Width, IDW, plus the following

measurement whenever applicable:

(i) 2" for fixed sides;

(II) 3" for hinged flip down sides;

(III) 4" for removable sides when sides not removed;

(IV) 4.5" for angle adjustable locking sides; and

(V) padding thickness

5. Outside Width is the Inside Width for removable sides when sides removed

6. Platform to Floor Height: measurement from the ground to the top of the platform below the ABS.

7. Maximum 4 sides/rails: head, foot, left side and right side.

8. Side/Rail Height: measurement from the top of sides to the top of steel platform, ABS, or base pads, whichever is the smallest

9. Factory setting of 8" Side/Rail Height is measured from top of steel platform to the top of sides

